



Federal TVET-System Medical Laboratory NTQF-III

Learning Guide:51

**Unit of competence: Performing
Urinalysis and body fluid analysis**

**Title: Performing Urinalysis and
body fluid analysis**


LG Code : HLT MLT3 M10Lo3-LG51

TTLM Code : HLT MLT3 M10 0919

LO3. Perform urinalysis tests





	Learning guide #52	Unit	Performing Urine and Body Fluid analysis
		Module	Perform urinalysis tests

Welcome to the module “Performing Urine and Body Fluid analysis”. This learner’s guide was prepared to help you achieve the required competence in “**Medical laboratory services Level-III**” this will be the source of information for you to acquire knowledge and skills in this particular occupation with minimum supervision or help from your trainer.

Summary of Learning Outcomes

After completing this learning guide, you should be able to:

- 3.1. Assembling required **equipment ,materials and systems**
- 3.2. Selection of the authorized tests
- 3.3. conduct Individual tests according to standards
 - 3.3.1 Physical Examination of Urine
 - 3.3.2 Chemical Examination of Urine
 - 3.3.3 Microscopic Examination of Urine
 - 3.3.4 Body Fluid Analysis
 - 3.3.4.1 CSF Analysis
 - 3.3.4.2 Semen Analysis
 - 3.3.5 Applying required quality control procedures
- 3.4 Recording interpretation of results
- 3.5 discussing of Colleagues with when result interpretation is outside parameters
- 3.6 verifying of Results before releasing for clinician/client
- 3.7 storage of Tested Samples and sample components for retesting when requested

Learning-instructions

1. Read the contents of this Learning Guide. It is divided into sections that cover all the skills and knowledge that you need.
2. Read the information written in the “Information Sheet #1, #2, #3, #4, #5, #6, #7, and # 8”.
3. Accomplish the “Self-check #1 on page 6, #2 on page 10, #3 on page 19, #4 on page 30, #5 on page 48, #6 on age 56, #7 on page



61

4. If you earned a satisfactory evaluation on self-check proceed to next learning Guide. However, if your rating is unsatisfactory, see your teacher for further instructions.
5. Read the “Operation Sheet” on page #31, #50, and #57, and try to understand the procedures discussed.
6. Practice the steps or procedures as illustrated in the operation sheet. Go to your teacher if you need clarification or you want answers to your questions or you need assistance in understanding a particular step or procedures

**Instruction Sheet #1****Learning Guide #49**

This learning guide is developed to provide you the necessary information regarding the

Following content coverage and topics –


LO1. Identify concepts of urinalysis

- 3.1. Assembling required **equipment ,materials and systems**
- 3.2. Selection of the authorized tests
- 3.3. conduct Individual tests according to standards
 - 3.3.1. Physical Examination of Urine
 - 3.3.2. Chemical Examination of Urine
 - 3.3.3. Microscopic Examination of Urine
 - 3.3.4. Body Fluid Analysis
 - 3.3.4.1. CSF Analysis
 - 3.3.4.2. Semen Analysis
 - 3.3.5. Applying required quality control procedures
- 3.4. Recording interpretation of results
- 3.5. discussing of Colleagues with when result interpretation is outside parameters
- 3.6. verifying of Results before releasing for clinician/client
- 3.7. storage of Tested Samples and sample components for retesting when requested

Learning Activities

1. Read the information written in the “Information Sheets”.
2. If you earned a satisfactory evaluation proceed to next module. However, if your rating is unsatisfactory, see your teacher for further instructions.
3. Read the “Operation Sheet” and try to understand the procedures discussed.
4. Practice the steps or procedures as illustrated in the operation sheet. Go to your teacher if you need clarification or you want answers to your questions or you need assistance in understanding a particular step or procedure



 Information sheet #1	LG #51	Perform urinalysis tests
	Topic	Assembling required equipment ,materials and systems

LO3. Perform urinalysis tests

3.1. Assembling required **equipment ,materials and systems**

1. The necessary materials used for the collection, centrifugation and examination of urine specimens are:

- ✓ Clean dry plastic or Glass containers, which enable to collect at least up to 15 ml of urine for routine urinalysis.
- ✓ Hand (manual) or electrical centrifuge.
- ✓ Conical centrifuge tubes, or regular test tubes.
- ✓ Pasteur pipette with rubber fit or automatic pipettes if possible.
- ✓ Slides and cover slides 20 x 20 mm.
- ✓ Electrical or solar microscope, which has 10x and 40x objectives.

2. Preparation of patient

- Explain the purpose of the test by using simple language. Do not use medical terms or try to explain details of the procedure.
- Advise the patient how to collect the specimen. The first morning urine or mid-stream urine specimen is more preferable, because it is more concentrated.
- If the patient is female, advise her to wash her genital organ before giving the specimen. This is because bacteria that are normally found on the genital tract may contaminate the sample and affect the result.
- Advise the patient to collect at least 15 ml of urine in to the clean, sterilize and dry urine cup that is supplied from the laboratory.

5.1.2. Source of Errors in the Microscopic Examination of Urine

Possible errors that may encounter during microscopical examination of urine include:

- Drying of the specimen on the slide.



- During trial of observing 2 specimens in a single slide by putting at each side of slide, (mix up of the specimens).
- If the supernatant fluid after centrifugation is not poured off properly, that is if some drop is left in the tube, it may decrease concentration of urine sediments and false result may be reported
- If the whole sediment with supernatant is discarded during inverting down the tube for long period, the whole sediments will be discarded and so again false negative result will be reported.



Self-check #1

Written tests

Instruction 1:- Say true or false for each of the following question

1. Explaining the purpose of the test by using simple language is better than using medical terms or try to explain details of the procedure.
2. If the whole sediment with supernatant is discarded during inverting down the tube for long period, may cause false negative result to be reported.
3. Drying of the specimen on the slide is the possible sources of error.

Instruction 2:- Answer the following question appropriately

1. List at least 3 sources of error during urine analysis?
2. List at least 3 materials used for examination of urine?

Answer sheet

Instruction 1

1. _____
2. _____
3. _____

Instruction 2

4. _____
5. _____

Score _____

Remark _____

Name _____ ID.no _____ Date _____



refractive index similar to that of urine in which they are suspended. Phase contrast has the advantage of hardening the outlines even the most ephemeral formed elements.

B. Stained Preparation

Cellular detail is best seen with stained preparation

The following stains are commonly used:

1. A crystal violet safranin stain (sternheimer and malbin) is useful in the identification of cellular elements.

Procedure

Add 1 or 2 drops of crystal violet safranin stain to approximately 1 ml of concentrated urine sediment. Mix and place a drop of this suspension on a slide and cover with cover slide.

Staining reaction to crystal – violet safranin stain:

RBC – Purple to dark purple.

WBC – Cytoplasm -violet to blue.

Nucleus – reddish purple.

Glitter cells – blue.

Cells	Nucleus	Cytoplasm
Squameous epithelial cells	Purple	Pink to violet
Euro epithelial	Dark blue	Blue
Renal tubular cells	Dark purple	Orange purple

2. Methyl blue (Loffler’s stain)

3. CytoDiachrome stains

When such stains are used, it is recommended that both the stained and unstained sediment be mounted and observed, as the stain may cause precipitation of some constituents. This is especially the problem with alkaline urine specimens, because the precipitated materials may obscure important pathological constituents.

II. Table 3. Relationship between Physiochemical and Microscopic Findings of Urine in Selected Disease States.

Physical Findings	Chemical Finding	Microscopic	Observation
Colored brown Turbidity Specific gravity	Protein + Blood +	WBC, RBC Hyaline or Granular or Cellular casts	Acute Glomerulonephritis
Urine volume Turbidity Odor pH	Protein + Blood + Nitrite +	- RBC, WBC - Cellular casts, - Bacteria	Acute tubular Necrosis Or lower Nephrosis
Specific gravity	Protein	Colorless	Cystinosis



Urine volume Specific gravity	Blood Protein + Blood +	Hexagonal Plate crystals	
Specific gravity Odor- sweet	Protein Glucose Ketone	Yeasts Some times Present	Diabetes Mellitus
Color darker	Glucose + Ketone + Blood + Bilirubin + Urobilinogen+	Pigment laden Prussian blue Casts	Hemochromatosis
Turbidity	Portion + Blood +	Casts Oval fat Bodies	Nephrotic Syndrome
Specific gravity	Protein + Blood +	Sickled RBC	SickleCell Syndrome
Turbidity	Protein +	RBC, WBC Casts	Systemic lupus Erythematosus

III. Table 4: Correct and Incorrect Approach in Urine Testing

Correct approach	Incorrect approach
Use fresh urine	Delay in the testing of urine without Preservation
Make quality control of reagents	Using expired reagents
Be aware of normal as well as abnormal results which are significant	Believing urine results have little significance in the overall diagnostic picture of the patient
Follow the directions carefully	Being careless
Accept only clear and proper collection bottles	Using any container.
Be familiar with interfering Substances	Not giving due attention to cross reaction and artifacts
Mix Urine properly	Not mixing well
Record results accurately	Not checking the results recorded during the training of new personnel
Give proper training to Professionals	New personnel always jumping into urinalysis because it is the easiest to do and least significant

**Self-check #2****Written tests**

Instruction 1:- choose the best possible answer for each of the following questions?

- Which of the following is the possible source of error in the microscopic examination of urine?
A. Inadequate centrifugation C. Inadequate centrifugation
B. Drying of the specimen on the slide D. All
- Identify the factors that may not results in falsely increase in high number of RBCs?
A. Menstrual bleeding C. Renal stone
B. Vaginal bleeding D. Aspirin ingestion or over dose
- Which one is not included in the chemical examination of urine?
A. Ketones B. Leukocyte Esterase C. Crystals D. Nitrite
- Creatinine clearance is most often used to monitor GFR because of the following reasons except?
A. Creatinine is an endogenous substance C. Creatinine freely filtered by glomerulus
B. Creatinine is reabsorbed by tubules D. None
- The correct method for labeling urine specimen containers is to.
A. Attach the label to the lid C. Attach the label to the container
B. Attach the label to the bottom D. Use only a wax pencil for labeling
- A urine specimen for routine urinalysis would be rejected by the laboratory because:
A. The specimen had been refrigerated C. The label was placed on the side of the container
B. More than 50 ml was in the container D. The specimen and its requisition did not match
- A sagittal/inner section of the kidney reveals three distinct regions. Which alternatives show those regions from outer to inner in sequences?
A. Pelvis→Cortex→Medulla C. Pelvis→ Medulla → Cortex
B. Cortex→Medulla→ Pelvis D. Medulla→ Pelvis→ Cortex
- Which of the following specimen type is/are commonly used for microbiological tests?
A. Random urine specimen C. Midstream urine specimen
B. Terminal urine specimen D. 24-Hour urine specimen
- A sagittal/inner section of the kidney reveals three distinct regions. Which alternatives show those regions from outer to inner in sequences?
C. Pelvis→Cortex→Medulla C. Pelvis→ Medulla → Cortex
D. Cortex→Medulla→ Pelvis D. Medulla→ Pelvis→ Cortex



10. The type of urine specimen that taken at any time of the day that the pts attend the diagnostic laboratory is termed as-----?

- A. Early morning specimen C. Random urine specimen
B. Midstream urine specimen D. Clean catch urine specimen

Note:- Satisfactory point is above five (>5)

Not- Satisfactory point is below five (<5)

Answer sheet

Instruction 1


1. _____
2. _____
3. _____
4. _____
5. _____
6. _____
7. _____
8. _____
9. _____
10. _____

Score _____

Remark _____

Name _____ ID.no _____ Date _____



 Information sheet #3	Unit	Performing Urine and Body Fluid analysis
	Topic	Physical Examination of Urine

3.1. Physical Examination of Urine

- Physical examination of urine is the first part of routine urinalysis.
- It is the simplest procedure of all urine examination, but this simplicity does not mean that anyone can do it without any background knowledge and experience.
- Physical examination of urine usually gives hint for the subsequent urinalysis.
- For example, white turbid urine sample may suggest to the technician the presence of Leukocytes (pus cells) and/or
- Epithelial cells in microscopic examination, and in chemical examination, with positive result of Nitrite.

3.1.1. Volume

Normally, 600 – 2000 ml of urine is voided per 24 hr.

Volume of urine excreted is related to:

- ✓ Individual fluid intake
- ✓ Body temperature
- ✓ Climate
- ✓ Individual's health status

Abnormally higher amount (greater than 2000 ml/24) or very low amount i.e. less than 600 ml/24hr occur mostly due to some pathological conditions.

For the measurement of the volume of urine, the patient should collect 24 hr urine specimen.

Clinical Significance

The Measurement of the volume of urine indicates the evaluation of fluid balance and kidney function.

When an individual excretes more than 2000 ml of urine/24 hr, consistently (for long period) it is called **Polyuria**.

It may occur due to:

- ✓ Diabetic mellitus
- ✓ Diabetic insipidus
- ✓ Certain tumors of brain and spinal cord
- ✓ Acromegaly
- ✓ Myxedema
- ✓ Some type of tubular necrosis (improper function of urine tubules)

Diuresis: Any increased amount of urine volume, even if for short period. It is usually due to excessive fluid intake.

Oliguria: Excretion of constantly small amount of urine, i.e. below 400 ml of urine/24 hr.

It may occur due to:

- ✓ Dehydration or poor blood supply to kidney that may be due to prolonged vomiting, diarrhea, etc.
- ✓ Obstruction of some area of the urinary tract/system (mechanical)
- ✓ Cardiac insufficiency



- ✓ Various renal diseases such as glomerulonephritis, etc.
- ✓ Fasting
- ✓ Excessive salt intake etc.

Anuria :Complete absence of urine excretion. It is less than 100 ml of urine per 24 hr. It may occur due to:

- ✓ Complete urinary tract obstruction
- ✓ Acute renal failure
- ✓ Acute glomerulonephritis
- ✓ Hemolytic transfusion reaction, etc

Polyuria: may result physiologically after consumption of

- ✓ Intravenous glucose or saline
- ✓ Coffee, alcohol, tea, caffeine
- ✓ Pharmacological agent, such as thiazides and other diuretics

2. Odor

Normally fresh voided urine from healthy individuals has faint aromatic odor, which comes from volatile acids, normally found in urine, mostly, ammonia.

The test is conducted by smelling of urine and the result is based on the perception of the technician.

Clinical Significance

Abnormal urine odor may result from aging of urine, disease and diet.

If the urine specimen is old, i.e. after collection, left on the bench without preservative for more than 2 hrs, it will have ammonical (pungent) odor.

The ammonical odor result is due to break down and conversion of urea in the urine into ammonia by the action of bacteria.

Cystinuria and homocystinuria (type of amino acids, voided from abnormal metabolism) have sulfurous odor.

Oasthouse urine disease has a smell associated with the smell of a brewery (yeast).

Tyrosenemia is characterized by cabbage like or “**fishy**” urine odor.

The presence of ketone bodies in the urine, that may be due to diabetes mellitus, vomiting, starvation, strenuous exercise, characterized by “**sweet fruity**” odor.

Butyric / hexanoic acidemia produce a urine odor resembling that of sweat.

Urine of infants, which has inherited amino acid metabolism disorder, smells like “**burnt sugar**” or maple, hence the name, “**maple sugar urine disease**”.

Also due to some food stuff such as asparagus, characteristic, urine odor is produced, which has no clinical significance.

3. Foam

Normally when urine specimen is voided in a container, it produces small amount of white foam. But during certain abnormal physiological and metabolic conditions, the color and amount of foam may be changed.

- ✓ For example, when there is high bile pigment in the urine, the amount of foam increases, and the color of foam becomes yellowish. This may indicate the presence of bilirubin in the urine.
- ✓ But the presence of yellowish foam should not be taken as a confirmatory test for the presence of bilirubin in urine. Chemical analysis of urine for bilirubin should be done.



4. Color

Normally color of urine may vary within a day; in the morning it has dark yellow color, while in the afternoon or evening, the color ranges from light yellow to colorless.

Normal urine color varies from straw (light yellow color) to dark amber (dark yellow).

- ✓ Light yellow indicate that the urine is more diluted, and has low specific gravity. Such exceptional condition occurs in case of diabetic mellitus. In this condition the color of urine is mostly light yellow, but because of having high glucose content, its specific gravity is high.
- ✓ On the other hand, dark amber (dark yellow) color mostly indicates that the urine is concentrated, and has high specific gravity. This type of urine is seen normally in the first morning urination.
- ✓ Normal urine color results from three pigments. They are:

1. **Urochrome**, responsible for yellow color formation. This pigments found in high proportion than the other two.

2. **Uroerythrin**, – responsible for red color formation.

3. **Urobilin**, – responsible for the orange-yellow color formation.

Thus, normal urine gets its color from a combination of the above-mentioned three pigments.

Procedure of the Test

Urine color is recorded, after looking at freshly voided urine specimen. If the urine sample color is not recorded within 30 minutes after collection, chemical changes will occur in it, and so its color will change, and will result in false report.

Clinical Implication

By observing the color of freshly voided urine, an experienced laboratory technician can forecast the possible findings in the chemical and microscopical examination of urine.

Depending up on the constituents of urine, the abnormal color of urine varies as follows:

- Pale to colorless urine may indicate:
 - Large fluid intake
 - Diabetic mellitus
 - Diabetic insipidus
 - Alcohol consumption
 - Nervousness
- Dark yellow or brown red urine may indicate:
 - Concentrated urine
 - Decreased fluid consumption
 - Dehydration
 - Fever
 - Certain urinary tract medication (e.g. phenazopyridine)
- Yellow brown or “beer brown” color may indicate the presence of bilirubin.

This is also confirmed:

- By looking at the yellow foam or green foam by shaking the sample.
- By letting it to stand for more than 30 minutes and looking at the change of color into green, because of oxidation of bilirubin into biliverdin.



- Due to bilirubin crystals, as mentioned in urine segment, the urine samples have opalescent appearance.
 - By doing chemical tests for bilirubin.
 - Clear red may indicate presence of Hemoglobinuria (presence of hemoglobin in the urine). This hemoglobinuria may result from:
 - Incompatible blood transfusion.
 - Increased red blood cell destruction (intravascular haemolysis) due to different hemoparasites, e.g. Malaria.
 - Glucose – 6-phosphate dehydrogenase deficiency.
 - Certain infections or disease.
 - Cloudy red / smoky red color may indicate hematuria (presence of red blood cell in the urine). It differs from clear red by the presence of RBC rather than Hgb alone. It is important to differentiate hemoglobinuria from hematuria, because the cause of this abnormal urine differs. On standing the red cells in hematuria may hemolyze and settle, and so the urine becomes clear red (hemoglobin in urine).
To differentiate this definition; specific gravity is important.
 - ✓ Hematuria has high specific gravity than hemoglobinuria.
 - Dark brown colored urine** may contain porphyrines, melanin, homogentistic acid, which is associated with an abnormal metabolism of tyrosine. Milky urine may contain fat, cystine crystals, and many WBC or amorphous phosphates.
- Dark reddish color** may indicate myoglobin (muscle Hgb), usually associated with extensive muscle injury, hemoglobinuria and porphyrine.

Interfering Factors

It is usually important to consider, that on standing of urine for more than 30 minutes, the urobilinogen that is found in urine will oxidize and change to urobilin. Thus due to this process, the color of urine becomes dark. Therefore, the physical examination of urine should be done immediately after the delivery of urine to the laboratory.

Other interfering factors that result in abnormal urine color formation are certain foodstuff, and medications.

- ✓ Food stuff, such as beets will give white red color.
- ✓ Drugs such as Vitamin B12 and riboflavin will give bright yellow color to urine.
- ✓ Rifampicin will give red color to urine.
- ✓ Iron salt will give dark color to urine.
- ✓ Sulfonamides will give rusty yellow or brownish color.

Therefore, when abnormal colored urine is observed, it is important to ask the patient, what kind of food he consumed in the last 36-24 hrs, and also whether he used drugs or not. If so, it is important to know what food and what drug he used.

5. Appearance (Transparency)

Fresh voided urine specimen is normally clear and transparent. On longstanding, due to chemical changes that occur in normal constituents of urine through time, it becomes turbid.

Procedure of the Test

- Appearance (transparency) of urine can be measured only by observation of fresh voided urine specimen.



□ Degree of cloudiness of the urine is described by using common terms, starting by clear to turbid i.e. clear, hazy, cloudy, very cloudy and turbid.

Clinical Implications

Freshly voided urine specimen appearance may indicate the presence of some abnormal constituents in it. Causes of turbid urine, as it is freshly voided include:

- White blood cells (pus cells) that occur due to UTI
- Kidney stones
- RBC's
- Yeast cells,
- High number of bacteria cells
- High number of epithelial cells
- Fat droplets in urine, which give opalescent appearance (rare condition).
- Amorphous urates, in case of gout and leukemia.
- High number of mucus trades.

All the above findings are confirmed by urine microscopic examination.

Interfering Factors

High consumption of foodstuff that contains urates and phosphates may produce cloudy urine. This is because of the precipitation of urates and phosphates in the form of amorphous urate and phosphates respectively.

Semen, or vaginal discharge mixed with urine is other common causes of urine turbidity. Urine specimen, stood for long period in the bench, will become hazy or cloudy due to precipitation of crystals, mucus trades etc., which normally occur in urine. The settlements of crystal and mucus trades seen in urine sample are to be preserved in refrigerator.

Amorphous urates have “**Brice red**” precipitation, while amorphous phosphates have **white** precipitations.

Clinical Significance

As indicated in the chapter one, one of the functions of renal system is to regulate pH of blood i.e. keeps pH of blood at 7.4 ± 0.05 . This is done by absorption or release of hydrogen ion, especially at distal convoluted tubules of the nephron, depending on the pH of blood, i.e. hydrogen ion absorbed from surrounding blood capillaries of nephron when pH is acidic (below 7.35), and release from nephron to the surrounding blood vessels when pH of blood is alkaline (above 7.45). pH measurement of urine, like other physical tests of urine, may indicate the on- going process in body, mostly about the renal system.

Normal pH of urine is 5-6.

***Persistent alkaline urine (pH > 6) may be caused by:**

- ✓ □ UTI
- ✓ □ Renal failure
- ✓ □ Vomiting
- ✓ □ Anorexia nervosa
- ✓ □ Alkalosis (metabolic or respiratory e.g. due to accumulation CO₂ in our body.
- ✓ □ Alkalizing drugs i.e. during intake of drugs such as streptomycin, kanamycin etc.



- ✓ It should also be important to bear in mind that certain vegetables, citrus fruits, and milk products also may cause alkaline urine, which is not pathological

*** Persistent acid urine (pH < 6) may be caused by:**

- Diarrhea
 - Malabsorption syndromes
 - Diabetic ketoacidosis
 - Dehydration
 - Fever
 - Starvation
 - And also certain drugs such as – Phenacetin
 - Here it is important to bear in mind that high protein diet may also result in acidic urine, but this is not a pathological condition.
 - pH measurement is also important in the management of renal stone patients, who are being treated for renal calculi and who are frequently given diets or medications to change the pH of the urine so that kidney stone will not form.
- Calcium phosphates, calcium carbonate, and magnesium phosphate stones develop in alkaline urine. In such instances the urine must be kept acidic (i.e. either by diet such as meat, or medication).
- Uric acid, cystine, and calcium oxalate stones are precipitated in acidic urine. Therefore, as part of treatment, the urine should be kept alkaline (either by diet e.g. leguminous plants, citrus fruits and most vegetables or by medication).

Interfering Factors

If urine specimen is left on the bench for more than 2 hours, bacteria will grow in it and by converting urea into ammonia, the pH will become alkaline. This is false alkaline urine, and indicates the specimen is not fresh urine.

6. Specific Gravity of Urine

Specific gravity is defined as the ratio of the weight of a fixed volume of solution to that of the same volume of water at a specified temperature, usually 20°C (in some books 25°C). The specific gravity of urine has been used for years as a measure of the total amount of material dissolved in it (total solids), and thus of the concentrating and excretory power of the kidneys.

Measurement of Specific Gravity

The following methods are used to test the specific gravity of urine:

- Urinometer
- Refractometer
- Reagent strip
- Weighing technique

Specimen: It should be the first urine passed at the beginning of the day with the patient having taken no fluid for 10 hours. The testing of random urine specimen has little clinical value.

1. The Urinometer

The specific gravity of a urine specimen is often measured with urinometer.

The urinometer is a glass float weighted with mercury, with an air bulb above the weight and a graduated stem on the top.



It is weighted to float at the 1.000 graduations in distilled water when placed in a glass urinometer cylinder or appropriate sized test tube. It is important that the cylinder, or test tube, be of the correct size so that the urinometer can float freely. The specific gravity of the urine is read directly from the graduated scale in the urinometer stem.

The scale of the urinometer is calibrated from 1.000-1.060 with each division being equal to 0.001.

Sources of Error:

- Temperature differences
- Proteinuria
- Glycosuria
- X-ray contrast media, it increases urine specific gravity
- Chemical preservatives

Urinometer Controls:

The following solutions can be used to check Urinometers:

Solutions Specific gravity pure water 1.000

Sodium chloride solution (2.5 g/dl) 1.018

“(5 g/dl) 1.035

“(7.5 g/dl) 1.051

2. Refractometer

It is an instrument, which reads the refractive index of the urine. There refractive index measurement depends on the number of dissolved particles in the urine.

The higher the concentration of the particles the greater the refractive index, and so the specific gravity.

3. Reagent Strip Test of the Specific Gravity of Urine

A test area to determine specific gravity in urine can be found in the multiple test strip of Ames called N-multistix. The reagent test area responds to the concentration of ions in the urine. It contains certain pretreated polyelectrolyte. The pKa of which changes depending up on the ionic concentration of the urine .The indicator bromothymol blue is used to detect the change. Colors ranges from deep blue when the urine is of low specific gravity through green to yellow- green when the urine is of high ionic concentration.



Self-Check 3	Written Test
---------------------	---------------------

Instruction1:- Say True or False

1. Urine color and urine concentration commonly vary together.
2. The normal yellow color of the urine is due primarily to uroblin, uroerythrin and urochrome.
3. A turbid urine specimen always indicates a pathologic condition.
4. The incidence of turbidity of the urine increases following refrigeration..
5. The pH of the urine usually rises after collection due to the growth of urea splitting bacteria, which produce ammonia.


Answer sheet

1. _____
2. _____
3. _____
4. _____
5. _____

Score_____
Remark_____

Name_____ID.no_____Date_____



 Information sheet #4	LG #51	Performing Urine and Body Fluid tests
	Topic	Chemical analysis of urine

3.3.2. Chemical analysis of urine

Chemical analysis of urine is an important procedure, in the detection of many diseases. Urine contains normal chemical compositions. But in abnormal conditions its composition varies in kind and quantities. So the chemical changes of urine can indicate disease at the early stage. The composition of urine varies because it is the principal route for soluble waste material from body metabolism. Its composition therefore depends greatly on how much and what specific waste material is to be excreted.

- Urea, creatinine, uric acid, ammonium salts, chlorides, sulphates and phosphates of sodium, potassium, calcium and magnesium are the **normal composition of urine**. They are excreted through the urine as a final body metabolism.
- Glucose, protein, ketone bodies, bilirubin, bile salts...etc are the **abnormal constituents of urine**. Normally these substances do not appear in the urine in detectable amount. So their appearance in the urine shows the pathological condition.
- For example, glucose does not appear in the urine in detectable amount. But during diabetes mellitus it appears in the urine. Protein also appears in the urine during renal disease. Generally the chemical examination of urine helps to investigate the health condition of individual.

1. Determination of Urinary Sugar (Glucose)

Glucose, a monosaccharide, is the principal sugar in blood, serving the tissues as a major metabolic fuel. It is mainly the end-product of carbohydrate digestion, which provides energy for life process. When body requires energy glucose oxidized to pyruvate and then to acetyl-CoA and enter cycle Krebs (tricarboxylic acid, TCA cycle). Along these metabolic processes it gives energy in the form of adenosinetriphosphate (ATP). ATP is very important energetic organic compound used for proper body function. When glucose is not required for the body's immediate energy needs, it is converted to glycogen and stored in liver and muscles by the metabolic process called glycogenesis. When there is an excess glucose in the blood (especially after carbohydrate meal), it can be also converted to fats. Glucose first oxidized to acetyl-CoA through glycolysis. The formed excess acetyl-CoA and then converted to fats to be stored in the tissue. When it is required to maintain the blood glucose level, particularly during starvation, glycogen is converted to glucose by glycogenolysis. For maintaining the blood glucose level, it can be synthesized from non-carbohydrate precursors like amino acids, glycerol, lactate and etc. by the metabolic process, which is called gluconeogenesis. The blood glucose level is controlled by a hormone, insulin, which is produced by the beta-islets of Langerhans of the pancreas. Insulin lowers the content of the glucose in the blood and increases its utilization and storage in the liver and muscle as glycogen. The absence or lower production of insulin resulted in Diabetes mellitus, which is characterized by an elevated blood glucose levels (hyperglycemia) and accompanying glycosuria and may be accompanied by changes in fat metabolism.



Glucose is the sugar most commonly found in the urine, although other sugars, such as lactose, fructose, galactose, and pentose, maybe found under certain condition. Normally, urine does not contain a sufficient amount of sugar to react with any of the popular enzyme or reducing tests. When sugar appears in the urine, it shows the abnormality caused by disease diabetes mellitus. Hence urine sugar tests are extremely useful in monitoring the treatment of diabetes.

Clinical Significance

The presence of detectable amount of glucose in the urine is known as glycosuria. Normally almost all the glucose, which passes from the blood into the glomerular filtrate, is reabsorbed back into the circulation by the kidney tubules (proximal convoluted tubules). Usually less than 15 - 20 mg/dl (0.8mmol) is excreted in the urine. But this amount cannot be detected by the routine laboratory tests. The term glycosuria is usually used to describe the presence of more than the normal amount (15- 20 mg /dl)of glucose in the urine.

The occurrence of glucose in the urine is not normal if more than 15 – 20mg/dl. The blood glucose concentration normally lies between 65 and 110 mg/dl. After a meal it may increase to 120 - 160 mg/dl. If the blood glucose concentration becomes too high (usually greater than 170 – 180mg/dl), the excess glucose will not be reabsorbed into the blood and glucose start appearing in urine. The lowest blood glucose concentration that will result in glycosuria is termed as the renal threshold. The most common condition in which the renal threshold for glucose exceeds is diabetes mellitus.

Causes of Glycosuria

- Physiological
- Pathological

1. Physiological

Sometimes under physiological situations, glycosuria can occur

- ✓ After large ingestion of carbohydrates
- ✓ Anything that stimulates sympathetic nervous system such as excitement, stress etc.
- ✓ 15 to 20% cases of pregnancy may be associated with physiological glycosuria.
- ✓ Renal Glycosuria: In some persons, glycosuria is found when blood glucose is in normal range. This is known as renal Glycosuria. This is again due to lowered renal threshold. Usually this is a benign condition.

2. Pathological Glycosuria

A. Diabetes mellitus

The most common condition for glycosuria is diabetes mellitus, a metabolic disorder due to deficiencies of insulin.

Glucose is not properly metabolized and blood glucose concentration rises, and when it is in range of 170 - 180 mg /dl , glucose starts appearing in urine.

B. Glycosuria due to other endocrine disorders

Deranged function of a number of endocrine disorders can cause hyperglycemia and this may result in glycosuria.

- e.g.
- Hyperthyroidism
 - Hyper adrenalism
 - Hyper pituitarism



Types of Urinary Sugar (Glucose) Tests

- Test for urine sugar is used to detect diabetes mellitus and also used to monitor the effectiveness of diabetic control.
- There are various tests for glucose which may be applied to urine.

The most frequently used are:

A. Non specific reduction tests based on the reduction of certain metal ions by glucose;

B. specific (Enzymatic) tests based on the action of glucose oxidase on glucose.

A. Non-Specific Tests for Glucose

These tests are based on the ability of glucose to act as reducing substances. Tests that are based on the reducing ability of glucose are not specific for glucose. In these tests, glucose is acting as a reducing agent, and any compound with a free aldehyde or ketone group will give the same reaction. Hence Glucose is not the only reducing substance that may be found in urine. Urine contains non-glucose reducing substances (NGRS) such as: uric acid, creatinine, galactose, fructose, lactose, pentose, ascorbic acid, chloroform, and formaldehyde.

Commonly used non-specific tests for urinary sugar are **Benedict's Qualitative Test** and the **Clinitest Tablet Test**.


1. Benedict's Qualitative Test

Benedict is a very sensitive copper reduction test and may give positive reactions with non-specific non-glucose reducing substances normally present in urine. Since glucose is the reducing agent, it is oxidized to gluconic acid. The positive reaction is indicated by a color change. It is a qualitative test in which the degree of color formation is proportional to the amount of reducing substance present in the specimen and the results are graded as negative, trace 1+, 2+, 3+, and 4+.

Principle

When boiled in an alkaline copper sulphate solution, glucose and other reducing substances reduce (convert) the blue copper (II) in Benedict's qualitative reagent to copper (I) oxide (Cu₂O), which is orange to red in color. A positive reaction is graded as a change in color ranging from blue to green, yellow, orange and finally red.

The overall reaction is:



The copper (II) ions are supplied in Benedict's qualitative reagent in the form of copper sulphate (CuSO₄). In the presence of a strong alkali this is converted to copper (I) oxide (Cu₂O). The heat is supplied by means of a boiling-water (100°C) bath. The tubes are brought back to room temperature, and the results are read when convenient.

Grade results according to the following criteria:

Negative: No change in the blue color of the reagent or the occurrence of a white or green precipitate from phosphates in the urine.

Trace: Slight amount of yellow precipitate with a greenish blue to bluish green mixed solution. (This represents less than 500mg/dl of sugar).



Principle

Diastix is a specific test for glucose based on the use of glucose oxidase, which is impregnated on the reagent strip. The chemical reaction is the same as for clinistix, the difference being the chromogen system used to indicate the presence of glucose. The reagent area contains glucose oxidase, peroxidase, a blue background dye, and potassium iodide as the chromogen. In a positive reaction oxidation of potassium iodide results in the formation of free iodide, which blends with the blue background dye to give shades of green through brown (The Boeringer dip-strip Test is also based on the same principle). As with clinistix, large amounts of ascorbic acid may give falsely negative or delayed results for glucose. This suppression is not as great as with clinistix, but it may cause problems.

Bleach and hydrogen peroxide may cause falsely positive reactions, as with Clinistix. Diastix has the advantage of being suitable as a screening test for the presence of glucose in the urine, and giving a rough estimate of the amount of glucose present. It detects as little as 100 mg of glucose per 100 ml of urine. However, urine specimens from pediatric patients must be subjected to a non-specific test for urinary sugar (Clinitest or Benedict's test) in addition to the specific glucose screening test in order to detect the presence of sugars other than glucose.

Sensitivity

Diastix reagent strip detects as little as 100mg of glucose in 100 ml of urine.

4.2 Determination of Ketone Bodies

Ketone bodies are normal products of fat metabolism. They are normally not detectable in the blood or urine. In normal metabolism, fat is broken down in the tissues to glycerol and fatty acids. The free fatty acids are transported by the plasma albumin to the liver where they are broken down to acetyl coenzyme A (acetyl Co-A) molecules. These condense with oxaloacetate in the Krebs cycle to produce citrate. The citrate is then oxidized to produce heat and energy. Whenever there is inadequate carbohydrate in the diet or a defect in carbohydrate metabolism or absorption, the body metabolizes increasing amounts of fatty acids, which is then converted into excessive amount of acetyl-CoA. The extra acetyl-CoA molecules join up in pairs to form acetoacetic acid. Most of this is reduced to β -hydroxybutric acid while some is decarboxylated to acetone. Acetoacetic and β -hydroxybutric acids are transported in the blood to the peripheral tissues to serve as an alternative fuel for cells. In the peripheral tissues these ketone bodies are reconverted to acetyl- CoA, and oxidized by the tri carboxylic acid cycle to give energy. Acetone is excreted in the urine.

Clinical Significance

When the rate of formation of ketone bodies is greater than the rate of their use, their levels begin to rise in the blood, which is called ketonemia, and eventually in the urine, which is known as ketonuria.

These two conditions are seen most often in cases of starvation and diabetes mellitus. Ketone bodies can be seen also in the urine during prolonged vomiting, severe diarrhea, anesthesia, severe liver damage, high fat intake and low carbohydrate diet. The excessive production and accumulation of ketone bodies may lead to ketosis.



Its physiological effect is serious because acetoacetic acid and β -hydroxybutyric acid contribute excess hydrogen ions to the blood, resulting in acidosis - a condition that tends to lower the blood pH. If not corrected in time this may result in death.

Another physiological effect of ketone accumulation concerns the substance acetone and acetoacetic acid. Both have been found to be toxic to brain tissue when present in increased amounts in the blood. So this condition can result in permanent brain damage.

When ketones accumulate in the blood and urine, they do not occur in equal concentrations. β -hydroxybutyric acid is present in the greatest concentration and acetone in the smallest concentrations. However most of the tests for ketonuria are most sensitive to the presence of acetoacetate. There are no simple laboratory tests for β -hydroxy-butricacid. Most tests react with acetone and acetoacetate or both.

Types of Tests for Ketone Bodies

A test for ketone bodies should be done routinely on any urine that is positive for glucose because they appear in the urine of diabetics. Test for ketones should be done within 2 hours after collection

Some of the commonly used tests for ketone bodies are the following:-

- Acetest tablet test,
- Acetone powder test,
- Reagent strip tests (Ex. Ketostix),
- Lang's test,
- Rothera's test.

Principle of the Tests

Both **acetone** and **acetoacetate** give a **purple color** with **alkaline sodium nitroprusside**. This is the general principle for the tests mentioned above.

Results - Report the test as positive or negative

4.3 Determination of Urinary Protein

Protein is a macromolecule, composed of one or more polypeptide chains, each possessing a characteristic amino acid sequence and molecular weight. It has many biologically important functions. Some of the functions are acting as enzyme(e.g. trypsin), transport protein (e.g. hemoglobin, myoglobin) nutrient and storage protein (e.g. ovalbumin(egg), casein (milk), contractile or motile protein (e.g. actin, myosin)structural protein (e.g keratin, fibroin, collagen), defense protein (e.g. antibodies, fibrinogen), and regulatory protein (e.g. insulin, growth hormone).

Test for urinary protein is one of the most important and valuable parts of the routine urinalysis. Albumin is one of the important proteins, which appears in urine during a pathological condition. It often occurs as a symptom of renal disease. Globulins are excreted less frequently. Bence Jones protein is a specific type of globulin excreted in multiple myeloma.



Clinical Significance

The presence of protein in the urine is called Proteinuria. It is one of the most important indicators of renal disease. Its presence in the urine depends on the nature of the clinical and pathological disorder and the severity of the specific disease.

Causes of Proteinuria

1. Increased permeability of the glomerulus

Normally, the glomerular membrane, the initial stage in the formation of urine, is not permeable for protein molecules. If the glomerular membrane is damaged these large protein molecules can pass through, and end up in the urine.

2. A decrease in normal re-absorption in the tubules

Under normal conditions, the small amount of protein (with lower molecular weight), which does filter through the glomerulus, is reabsorbed back into the blood stream. Normal urine, therefore, contains only traces of protein, insufficient for detection by routine laboratory tests. However, the concentration of protein that normally filters into the glomerular filtrate is extremely small, and only 1% of the glomerular filtrate is eliminated from the body as urine; the rest is reabsorbed. Failure to reabsorb any protein from this large volume of glomerular filtrate will result in fairly large amounts of protein in the urine.

Types of Proteinuria

1. Accidental or false proteinuria

Accidental or False Proteinuria occurs when there is a mixture of urine with a proteinous fluid such as pus, blood or vaginal discharge. These can occur in infection of the kidney, bladder or vagina.

2. Physiological or functional proteinuria.

Physiological or functional proteinuria is protein excretion in association with fever, exposure to heat or cold, excessive exercise, emotional stress, and later stage of pregnancy. The underlying physiologic mechanism that induces proteinuria in all of these, is renal vasoconstriction.

3. Postural (orthostatic) proteinuria

Postural or orthostatic proteinuria is excretion of protein by patients, who are standing or sitting for a longtime. The proteinuria is intermittent and disappears when the individual lies down. It can also occur during abnormal curvature of spinal cord.

4. Renal or true proteinuria

Renal or true proteinuria occurs when protein passes from the blood in to the urine because of some malfunction in the filtering system, either in the glomerulus or tubules.

Table .2 Proteins in Urine

Proteins	Conditions
Albumin	<ul style="list-style-type: none"> ✓ Strenuous Physical Exercise ✓ Emotional Stress ✓ Pregnancy ✓ Infections ✓ Glomerulonephritis ✓ Newborns (First Week)
Globulin	<ul style="list-style-type: none"> ✓ Glomerulonephritis



	✓ Tubular Dysfunction
Hemoglobin	✓ Hematuria ✓ Hemoglobinuria
Fibrinogen	✓ Severe renal disease
Nucleoprotein	✓ WBCs in Urine ✓ Epithelial Cells in Urine
Bence jones	✓ Multiple Myeloma ✓ Leukemia

Tests for Urinary Protein

A. *Precipitation or Turbidimetric Tests*

Principle: The general principle of these tests is that protein is either precipitated out of the urine specimen by means of a chemical, which is usually a strong acid, or it is coagulated out of solution with heat. These tests include:

- Robert's test
- Heller's test
- Sulphosalicylic Acid Test & Heat and Acetic Acid Test

Turbidimetric test based on acid reagents are non-specific since any urine components, which is insoluble in acid, will give a positive result.

It requires large volumes (0.5 to 5 ml) and requires either disposable tubes or glass tubes which must be cleaned for re-use.

The results of the precipitation tests are read in terms of the amount of precipitate or turbidity that is formed in a test tube (in case of Heat and acetic acid, and Sulphosalicylic acid tests) or in terms of the size of ring of contact between reagents in case of Robert`s and Heller`s tests. The amount of turbidity or precipitation is roughly proportional to the amount of protein present in the urine specimen, and the results are generally graded as negative, trace, 1+, 2+, 3+, or 4+.

Since the result in precipitation tests is determined by the presence of either turbidity or a precipitate, it is important that the urine be free from particles or clear before the test is performed. To clear the urine, it should be filtered or centrifuged. The clear filtrate is tested for the presence of protein.

The **non-ring** precipitation is read and interpreted as follows:

Negative - no turbidity or no increase in turbidity (approximately 5mg/dL or less)

Trace - Perceptible turbidity (approximately 20 mg /dL).

1+ - Distinct turbidity, but no discrete granulation (approximately 50 mg/dL).

2+ - Turbidity with granulation, but no flocculation (approximately 200 mg/dL).

3+ - Turbidity with granulation and flocculation (approximately 500 mg/dL).

4+ - Clumps of precipitated protein, or solid precipitate (approximately 1000mg/dL or >)

The Ring Test is read as follows:-

Negative - No cloudiness appears at the zone of contact

Trace - Ring is just perceptible against a black background

1+ - Ring is distinct against a black background, can barely be seen when held up to the light.

2+ - Ring is very definite against light, fairly visible when viewed from above



3+ - Ring is heavy against light, distinct cloudiness when viewed from above.

4+ - Ring is thick and dense against light, opaque when viewed from above.

The reading is interpreted as in the case of non-ring precipitation test.

A. Robert's Test

Principle

The principle of this test is based on the precipitation of protein and formation of white compact ring using concentrated Nitric acid(HNO₃).

C. Sulphosalicylic Acid Test

Principle

This test is based on the precipitation of protein (particularly albumin) by sulphosalicylic acid,

D. Heat and Acetic Acid Test

Principle

The test is based on the precipitation of protein by heat.

Sensitivity

This method is the most sensitive for small amount of protein and can reliably detect protein concentrations of 2 to 3 mg/dl.

II. Colorimetric Reagent Strip (Dipstick) Tests

The Colorimetric (dipstick) Protein Tests are more specific than Turbid metric Tests.

They require only a drop of urine enough to moisten the reagent area. The Colorimetric reagent strip test is based on the ability of protein to alter the color of some acid-base indicators without altering the pH. When an indicator, such as tetra bromophenol blue is buffered at pH 3, it is yellow in solutions without protein but , in the presence of protein, the color will change to green and then blue with increasing protein concentrations. In this case the pH of the urine is held constant by means of a buffer so that any change of color of the indicator will indicate the presence of protein.

The tests for urinary protein are all commercial ones that are available as reagent strip, tests (Dipsticks) either alone or in combination with other tests. Example. Albustix, Uristix, N-Multistix, Combur3 or Combur9. Although the colorimetric tests are useful primarily as screening tests for protein, these strip tests can be read semi quantitatively as negative, trace, 1+, 2+, 3+, or 4+ to give a rough estimate of the amount of protein present. To do this, the resulting color must be matched closely with the color chart provided with the test strips. The albustix and other multiple-reagent strips produced by amesco. are plastic strips with protein test areas impregnated with citrate buffer and tetra bromophenol blue. The citrate buffer maintains the pH at 3. At pH 3 tetrabromophenol blue is yellow in the absence of protein and yellow - green, or blue in its presence. The shade of the color is dependent on the amount of protein present. Falsely positive reactions may occur when protein is absent, if the urine is exceptionally alkaline or highly buffered.

Quantitative 24 hour Protein Determinations

Simple estimates of the protein content of urine are performed by quantitating the amount of precipitation formed following the addition of a specific chemical to the urine. The precipitate is measured either by comparison with known standards (sulphosalicylic acid turbidity test) or by recording the height of the column of precipitate in a specially designed tube (Esbach's test).



4.4 Determination of Bilirubin

Bilirubin is a waste product that must be eliminated from the body. It is formed by the breakdown of hemoglobin in the reticulo-endothelial cells of the spleen and bone marrow, and then transported to the liver.

On its way to the liver it is not water-soluble, and is carried through the blood stream linked to plasma albumin. This water insoluble form of bilirubin is often referred to as free bilirubin or unconjugated bilirubin or indirect bilirubin. Since this albumin - bound form is insoluble in water; it does not appear in the urine.

In the liver bilirubin is converted to a water soluble product by conjugation with glucuronic acid to form bilirubin glucuronide. The water-soluble form is called conjugated bilirubin. It is also called direct bilirubin. The liver cells that form the conjugated bilirubin excrete it into the bile and it is then excreted into the intestinal tract through the bile duct. In the small intestine this conjugated bilirubin is converted by intestinal bacteria to urobilinogen or stercobilinogen. Even though normally the level of conjugated bilirubin in the blood is not high enough to cause significant amounts to appear in the urine, this water soluble and conjugated bilirubin can be excreted by the kidneys.

Normal Value: approximately up to 0.02 mg/dl (This amount is not detected by routine qualitative or semi quantitative techniques).

Clinical Significance

Tests for urinary bilirubin and urobilinogen were normally performed only indicated by abnormal color of the urine or when liver disease or a hemolytic condition was suspected from the patient's history. The presence of bilirubin and urobilinogen in the urine is an early sign of liver cell disease (hepatocellular disease) and obstruction of the bile flow from the liver (Obstructive or post - hepatic jaundice).

Urine containing bilirubin will typically have been brown color and produce a yellow foam when shaken. Bilirubin is not stable in solution, but will be oxidized to biliverdin, which is a green pigment. Thus urine containing bilirubin will typically be red-brown when voided, and will turn green on standing, especially if exposed to light. Tests for bilirubin will not be positive in the presence of biliverdin; so the urine must be examined when fresh.

Tests for Bilirubin

Tests for bilirubin are based on the oxidation of bilirubin to biliverdin.

Specimen: Freshly passed urine is required. Urine containing bilirubin should be analyzed immediately after collection (with in 2 hrs of voiding). If bilirubin exposed to sunlight, it will oxidize to biliverdin, which cannot be detected by the reagents used in any of the tests. The following tests are used to detect bilirubin in the urine.

C. Diazotization Tests for Bilirubin

The tablet and reagent strip tests for bilirubin are based on the coupling of bilirubin with a diazonium salt in an acid medium to form azobilirubin, which gives a blue or purple color.

1. Icotest Tablet Test

The Icotest tablet contains nitrobenzinediazonium, p-toluene sulfonate (bilazo), sulfosalicylic acid, and sodium bicarbonate. The mats are absorbent as best as cellulose.



2. Reagent Strip Tests for Bilirubin (Ex. Multistix)

Principle

These tests for bilirubin are available only on multiple-reagent strips in conjugation with other tests. They are diazotization tests and are analogous to the Ecotest tablet test. The test area for bilirubin on Multistix and other Ames Co. reagent strip products is impregnated with 2,4-dichloro-aniline diazonium salt. The reagent strip tests for bilirubin are difficult to read and the color formed after reaction with urine must be carefully compared with color chart supplied by the manufacturer.



Self-Check 4	Written Test
---------------------	---------------------

Answer the following questions

1. Discuss by comparison the Benedict's Qualitative and Glucose oxidase Tests.
2. List down the possible substances, which give false positive results in non-specific tests for glucose determination.
3. Mention the physiological effects of ketone accumulation in blood.
4. Write the principle of the test for determination of bilirubin and hemoglobin.
5. Write the general principles for the two types of determination of urinary protein.


Answer sheet

1. _____
2. _____
3. _____
4. _____
5. _____

Score _____
Rate _____
Remark _____

Name _____ Date _____



 Operation sheet #2	LG #51	Performing Urine and Body Fluid analysis
	Topic	Performing Chemical Examination of Urine

4.1. Procedure of performing chemical urinalysis by reagent strips

Observe the precautions and follow the instructions supplied by the manufacturer.

1. Wear gown, glove and other PPE
2. Clean the working bench
3. Assemble the required materials
4. Collect 15ml of urine into clean, dry container
5. Dip the reagent area of the strip briefly into the specimen.
6. Remove excess urine by tapping or drawing the edge of the strip along the rim of the urine container.
7. Compare the color that develops with the color chart supplied by the manufacturer and report as indicated on the chart.

4.2. Quantitative 24 hour Protein Determinations

Purposes

Simple estimates of the protein content of urine are performed by quantitating the amount of precipitation formed following the addition of a specific chemical to the urine. The precipitate is measured either by comparison with known standards (sulphosalicylic acid turbidity test) or by recording the height of the column of precipitate in a specially designed tube (Esbach's test).

Procedure

- a. Pipette 2.5 ml of centrifuged urine into a test tube.
- b. Add 7.5 ml of 3% sulphosalicylic acid.
- c. Invert to mix
- d. Let stand 30 minutes.
- e. Compare the turbidity with known standards prepared from solutions containing 10, 20, 30, 40, 75 and 100mg albumin/dl, and estimate the concentration of the unknown. If the unknown urine contains more than 100mg/dl protein, dilute the urine and repeat the test.

LAP TEST	Practical Demonstration
Name_____	ID.No_____Date_____
Time started_____Time ended_____	


Instruction1:- Demonstrate each of the following activities.

Project1:- Performing urine chemical examination

Task1:- Perform urine dipstick tests of glucose determination?

Task2:- Perform Quantitative 24 hour Protein Determinations?



 Information sheet #5	LG #51	Performing Urine and Body Fluid analysis
	Topic	Microscopic Examination of Urine

5.0. Microscopic Examination of Urine

Microscopic examination of urine is one of the routine tests of urinalysis.

Urine contains many substances in addition to water.

The amounts of solid substances, which are found in the urine, may indicate an individual's health status. i.e. whether one is healthy or sick.

Normally small amount of solid substances is found in the urine. But when their concentration become high, it may indicate the existence of abnormal physiological function of our body. Microscopic examination of urine to some extent can be considered as **“renal biopsy”** because it reveals more about the function of the kidneys.

Repeated evaluation of urine sediment is frequently valuable in following the course and management of urinary tract disorders, because the appearance of cellular elements, and casts in the urine is a reflection of changes that take place in the kidney.

Urine sediments can grossly be categorized into organized and non-organized sediments based on the substances they are composed of.

Urinary Sediments

Classification of Urinary Sediments

Organized Elements

- Formed from Living Materials

Non-organized Elements

- Formed for Non-living Material (Crystals)

Organized (Formed) elements

- WBCs/HPF -Amorphous Urates,
- Epithelial cells / LPF -Uric acid crystals,
- Casts / LPF -Cystine crystals
- Parasites/LPF -Calcium Phosphate
- Bacteria / HPF -Cholesterol
- Ammonium Biurates
- Yeast Cells / LPF -Tyrosine, Leucien, Bilirubin,
- Mucus trade/LPF -Calcium sulfates (urates)
- Spermatozoa - Calcium carbonate
- Miscellaneous substances (Common contaminants)

Non-organized (Non-living Material)

I. Slightly acidic urine

- Triple phosphates
- Amorphous phosphate
- Calcium carbonate
- Calcium phosphate

II. Acidic, Neutral, or slightly Alkaline Urine crystals

- Calcium Oxalate crystals



III. Alkaline, Neutral, or Slightly acidic urine

- Triple phosphates

IV. Alkaline Urine Crystals

- Amorphous phosphate
- Calcium carbonate
- Calcium phosphate

5.4 Organized Urinary Sediments

A. RED BLOOD CELLS

Appearance: Normally RBCs appear in the fresh sample as intact, small and faint yellowish discs, darker at the edges

- Measure 7-8 μm
- In concentrated urine may be crenated, and their size became small (5-6 μm)
- In diluted urine, RBCs may be turgid and increase in size (9-10 μm)
- In alkaline urine, they may be small or entirely destroyed forming massive of brownish granules




Fig. 3.1. shows RBCs and Calcium oxalate

Clinical Implications: When the number of RBCs is found more than their normal range, usually greater than 5 RBCs/HPF it may indicate:

□ Presence of disease conditions in the urinary tract, such as:

- Acute and chronic glomerulonephritis
- Renal stone
- Cystitis
- Prostates
- Trauma of the kidney
- Presence of parasites, such as: schistosoma.
- Presence of bacterial infection, such as: renal tuberculosis
- Other disease conditions, such as hemophilia, malignant hypertension.

Temporarily (transient) increased RBC may be seen

- After strenuous exercise



- Exposure to cold temperature

Other substances confusing with RBCs

Yeast cells, and fat droplets may confuse with RBCs morphologically. They may be differentiated by their morphology.

Red blood cells are somewhat round or disc shaped, and uniform in size: while yeast cells are oval in shape, and have budding at the surface. On the other hand fat droplets are irregular in size and they are shiny.

Another means of differentiating RBCs from yeast and fat droplets is that, when 5% of acetic acid is added under the cover slide, RBCs will hemolyze, while yeast cell and fat droplets will not show any change.

How to report result:

- After looking RBCs under the 40x objective, they can be reported by mentioning the average number of RBCs/HPF.

Interfering factors:

Factors that may result falsely in high number of RBCs, i.e. without the presence of actual renal or other normal physiological disturbances included:

- Menstrual bleeding
- Vaginal bleeding
- Trauma to per anal area in female patients
- Following traumatic catheterization
- Due to some drugs, such as,
 - Aspirin ingestion or over dose
 - Anticoagulant therapy over dose

B. LEUKOCYTES (WBCs)

Normal range: 0-4 WBC/HPF.

Appearance: normally, clear granular disc shaped,

- Measure 10-15 μm , the nuclei may be visible.
- In alkaline urine, they may increase their size and become irregular.
- Predominantly, polymorph nuclear neutrophils are seen.
- Sometimes because of predominance of neutrophils and the occurrence of bacterial cell together with polymorph nuclear cells, WBCs are called pus cells.
- WBCs (pus cells) may be seen in clumps.
- It is also possible to see single irregular nuclei and small round lobed nuclei in the WBCs, that are seen in the urine sediment.




Fig. 3.2. WBCs

Clinical implication: Increased numbers of leukocyte urine are seen in case of:

- Urinary tract infection
- All renal disease
- Bladder tumor
- Cystitis
- Prostates
- Acute or chronic bacterial infection such as renal tuberculosis, temporarily increased numbers of leukocytes are also seen during: Fever and after strenuous exercise

How to report the result:

- After observing the distribution of leukocytes under 40 x objectives, at least 10 fields of microscope, it is possible to report as: 0-5leukocytes / HPF, 20-39 leukocytes / HPF or 0-5 leukocytes / HPF are seen..... normal

5-10 leukocytes / HPF are seen..... few leukocytes / HPF

10-20 leukocytes/HPF are seen.....moderate leukocytes/ HPF

20-30 leukocytes /HPF are seen many leukocytes / HPF

Above 30 leukocytes / HPF / are seen full/field

C. EPITHELIAL CELLS

• Normally few epithelial cells (0-2 / HPF) can be found

• Appearance

Their size differs depending on the site from which they originated.

a. renal cells

- Size is small as compared to other epithelial cells

- It measures 10 μ to 18 μ m in length, i.e., slightly larger than leukocytes

- Very granular

- Have refractive and clearly visible nucleus

- Usually seen in association with proteins or casts (in renal disease).

b. Cells from pelvis and urethra of the kidney

- Size is larger than renal epithelia's



- Those from pelvis area are granular with sort of tail, while those from urethra are oval in shape
- Most of the time urethral epithelia is seen with together of leukocytes and filaments (mucus trades and large in number)
- Pelvic epithelia's seen usually with no leukocyte and mucus trade, and are few in number

c. Bladder cells

- Are Squameous epithelial cells?
- Very large in size.
- Shape seems rectangular and often with irregular border.
- Have single nucleus.

* Here it is important to keep in mind that it is not expected from an experienced Lab. technician after simply observing epithelial cells, to say that these are urethral cells, and of pelvic origin and reporting such a false result in the laboratory request form.

* Knowing the origin of the epithelial cells and reporting it, may have more meaning when requested by the physician for special purpose, especially by the urologists.



Fig. 3.3 epithelial cells

Clinical implication

Presence of epithelial cells in large number, mostly renal types may indicate:

- Acute tubular damage
- Acute glomerulonephritis
- Silicate over dose

* The presence of large number of epithelial cells with large number of Leukocytes and mucus trades (filaments) may indicate Urinary tract Infections (UTI).

Reporting of the result:

- Epithelial cells distribution reported after looking under 10x (low power objective) of the microscope.
- Usually they are reported semi quantitatively by saying
 - Occasional epithelial cells /LPF 1-3 epithelial cells seen in the whole LPF
 - Few epithelial cells / LPF..... 2-4 epithelial / LPF
 - Moderate epithelial cells / LPF..... 6-14 epithelial / LPF
 - Many epithelial cells / LPF..... 15-25 epithelial/ LPF



- Full of epithelial cells / LPF.....when the whole field of 10 x objective covered by epithelial cells.

CASTS

• Formed by precipitation of proteins, and aggregation of cells within the renal tubules. Most of them dissociate in alkaline urine, and diluted urine (specific gravity ≤ 1.010) even in the presence of Proteinuria. Most of them are transparent. Thus to look them clearly, it is important to lower the condenser and close (partially) the diaphragm. Look them under 10 x (low power objective) of the microscope. There are different kinds of casts based on their shape and content (morphologically) may be grouped in to the following.

A. Hyaline Casts

• Normal range: 0-2/HPF

Appearance

- Transparent (clear), cylindrical shape
- Have parallels side with slightly round ends
- Their appearance in urine depends on rate of urine flow, i.e. many hyaline casts are seen when the flow rate is slow, and are not seen in alkaline urine mostly; and as the degree of protein urea is high, there concentration also increase.

Clinical Implication

Presence of large number of hyaline casts may show possible damage of glomerular capillary membrane. This damage permits leakage of protein through glomerulus and result in precipitate and gel formation(i.e. hyaline casts) in the tubule. Thus this may indicate:

- Nephritis
- Meningitis
- Chronic renal disease
- Congenital heart failure
- Diabetic nephropathy

Hyaline casts may also be seen in moderate number temporarily in thecase of:

- Fever
- Postural orthostatic strain
- Emotional stress
- Strenuous exercise
- After anesthesia

B. Granular Casts

• More similar in appearance with hyaline casts and in which homogenous, course granules are seen. More dense (opaque)than hyaline cast, thus can be more easily seen than hyaline casts. They are also shorter and broader than hyaline casts. May represent the first stage of epithelial cell cast degeneration.

Some other studies also suggest that, they are formed independently from cellular cast degeneration, and stated that they result from aggregation of serum proteins into cast matrix of mucoproteins

• Based on the amount and type of granules, they can be further divided into fine, and course granular casts.

Clinical implication



Granular casts may be seen in

- Acute tubular necrosis
- Advanced granulonephritis
- Pyelonephritis
- Malignant nephrosclerosis
- Chronic lead poisoning
- In healthy individuals these casts may be seen after strenuous exercise

C. Cellular Casts

Cellular casts are casts, which contain

- Epithelial cells
- White blood cells
- Red blood cells

Normal range: normally not seen in normal individual

Appearance

- These are casts in which cellular elements are seen.
- Formed usually after accumulation of cellular element in the renal tubules

Clinical Significance

- Epithelial / renal / casts mostly seen in tubular degeneration.
- Red cell cast usually seen in acute glomerulonephritis cases.
- White blood cell casts seen mostly during pyelonephritis conditions.

NOTE: Casts are very significant findings of urine microscopic examination. This is because their presence indicates the existence of renal disease. Sometimes it is possible to get a single cast having coarse granules, fine granules and fat droplets, i.e. different substances in a single cast, at the same time. At this time decision is made after looking and evaluation of other fields and based on the majorities.

Reporting of Laboratory Result

- Casts are examined under 10x objective of the microscope.
- Always the condenser should be lowered and at the same time in order to have good contrast, the diaphragm should be partially closed.
- Casts are reported quantitatively by saying:

- o Occasional casts / LPF
- o Few casts / LPF
- o Moderate casts / LPF and
- o Many casts / LPF

During reporting the type of cast that is seen should also be mentioned

Example: few hyaline casts / LPF are seen

PARASITES

Parasites that can be seen in urine microscopy are:

- Trichomonas vaginalis
- Schistosoma haematobium
- Wuchereria bancrofti

* Other parasites also may occur due to contamination of the urine with stool.

A. Trichomonas Vaginalis

It is a protozoa parasite that infects the genitourinary tract.

Appearance

Medical Laboratory Level III	Vision :01 Sep. 2019:	Page 41 of 70
	Copyright Info/Author: Federal TVET Agency	



- Size is about 15 μm .
- Shape is round, globular.
- Has vibratory, whirls and turns type of movement.
- Has also undulating membrane that is like the fin of a fish, on one side very motile.
- Have 4 flagella.




Figure 3.4. Trophozoites of *T. vaginalis*

B. Schistosoma Haematobium

It is fluke that infect venules of the bladder.

Appearance of the egg

- It is found in the urine sediment.
- Has pale yellow brown color.
- Large and oval in shape.
- Has characteristic small spine at one end (terminal spine).
- Measure about 145 x 55 μm .
- The egg contains a full-developed miracidium. Sometimes the miracidium hatch from the egg and can be seen swimming in the urine. The miracidium swim in the urine by the help of ciliates that are surrounding it. High excretion of *S. haematobium* egg can be seen usual between 10.00a.m. and 2 p.m. It is also important to remember that even when persons are highly infected, eggs may not be present in the urine. Therefore that is important to examine several specimens collected on different days and examine carefully, that is due to the irregular pattern of egg excretion.



Figure 3.5. Egg of ***Schistosoma Haematobium***



C. Wuchereria Bancroftie

- It is tissue nematode that invades lymph vessels. It is usually attack lower limb.
- In chronic bancroftie filariasis, a condition called chyluria can occur. i.e. passing of chyle in the urine. It occurs when the urogenital lymphatic vessels, which are linked to those, that transport chyle from the intestine became blocked and rupture.
- Chyle consists of lymph and particles of digested fat (soluble in ether).
- Urine containing chyle appears creamy white. When blood is also present, the urine appears pinkish-white.
- Large, measuring 275-399 x 8-10 μm .
- Body curves are few, nuclei are distinct.
- Sheath stains pink with Giemsa and palely with haematoxylin.
- There is no nuclei in the tip of at the tail.

Other points that should be considered also

- The parasite usually found in high concentration during night from 10:00 p.m. – 4:00 a.m. and i.e. it has nocturnal periodicity.
- Differentiate from *B. malai* and *L. loa* by its tail feature.
- Differentiate from *Mansonella* species by its large size and sheath.

YEAST CELL

Yeast cells are fungi that are not normally seen in health individuals.

Appearance

- Variable in size
- Colorless.
- Oval in shape, and usually form budding.
- Have high refractive index.
- Usually confused with Red Blood Cells. The way in which one can differentiate yeast cells from RBC is discussed in detail under Red Blood Cells.




Figure 3.6. budding yeast




Figure 3.7. branching pseudohyphae

Clinical Significance

- They are usually of *Candida* species (*Candida albicans*) and are common in patients with
 - Urinary tract infection
 - Vaginites
 - Diabetic mellitus
 - Intensive antibiotic or immunosuppressive therapy.

BACTERIA

Bacteria are the most common cause of UTI and aerobic gram-negative bacilli, particularly, members of the enterobacteriaceae, are the most dominant agents. The Gram-positives account for proportionately large number of infections in hospital inpatients. Normally, bacteria are not seen in the healthy individual's urine.

To check the presence or absence of bacteria a technician can either check for Nitrate that was formed in the urine after breakdown of nitrite into nitrate by the metabolic action of bacteria. Hence, dipstick test can give indirect clue. Or one can use urine microscopy test to check the presence of pus cells within the drop of urine or its sediment. Further the observed bacterial cell can be identified by bacteriological culture.

Appearance

- Bacteria that are seen in the microscopic examination of the drop of urine sample. Their shape varies with the type of bacteria observed..
- Depending on the type of bacteria they can be either motile or non motile organisms.
- They can be observed when examined under less than 40 x (high power) objective of the microscope.

Clinical Significance

- Presence of bacteria may indicate the presence of UTI or contamination by genital or intestinal micro flora.
- To confirm what type of bacteria they are and whether or not they are the causes of the disease, it is important to culture them in appropriate media and perform biochemical tests for identification.

Report of the Result

The bacteria concentration before or without performing culture and identification of the bacteria can be reported as:



- Occasional bacteria / HPF
- Few bacteria / HPF
- Moderate bacteria / HPF
- Many bacteria / HPF
- Full of bacteria / HPF.

Elements in Urinary Sediment	Usual Distinguishing Color of Stained Elements	Comments
Squamous epithelial cells	Dark shade of orange-purple	Light purple or blue
Inclusions and Matrix		
Hyaline casts	Pale pink or pale purple	Very uniform color; slightly darker than mucous threads
Coarse granular inclusion casts	Dark purple granules in purple matrix	
Finely granular inclusion casts	Fine dark purple granules in pale pink or pale purple matrix	
Waxy casts	Pale pink or pale purple	Darker than hyaline casts, but of a pale even color; distinct broken ends
Fat inclusion casts	Fat globules unstained in a pink matrix	Rare; presence is confirmed if examination under polarized light indicates double refraction
Red cell inclusion casts	Pink to orange-red	Intact cells can be seen in matrix
Blood (hemoglobin) casts	Orange-red	No intact cells
Bacteria	Motile: do not stain Nonmotile: stain purple	Motile organisms are not impaired
<i>Trichomonas vaginalis</i>	Light blue-green	Motility is unimpaired in fresh specimens when recommended volumes of stain are used; immobile organisms also identifiable
Mucus	Pale pink or pale blue	
Background	Pale pink or pale purple	

5.5 Non-organized Elements (Urine Crystals)

Appear usually after the specimen (urine) collected and left without examination. Mostly occur during metabolic abnormalities and excessive consumption of certain foodstuffs. May be classified into acidic, basic, and both acidic and basic based on:

- PH of urine in which they are usually seen.
- Solubility characters.

Identification of particular urine crystals from patient urine-sediment mainly serves as

- Guide to diagnose most likely type of calculus present.
- Mode of therapy of calculus by adjusting of urine, and by avoiding the intake of certain calculus precursors.



- Occurrence of certain abnormal urine crystals, such as cystine, Leucine, and Tyrosine, indicate the patient is in certain metabolic disorders and some drug crystals in the urine include sulfonamides, aspirin, and caffeine, used to follow the treatment condition.

The Usual Crystals Found in Urine	
Alkaline pH	Acid pH
Amorphous phosphates	Amorphous urates
Triple phosphates	Uric acid
Ammonium biurates	Calcium oxalates
Calcium phosphates	
Calcium carbonates	Cystine

Normal Crystals

- Uric acid Crystals
- Calcium Oxalate Crystals
- Hippuric Crystals
- Calcium Phosphate Crystals
- Triple Phosphate Crystals
- Calcium Carbonate Crystals
- Ammonium Biurate Crystals

Abnormal Crystals

- Bilirubin Crystals
- Cholesterol Crystals
- Cysteine Crystals
- Leucine Crystals
- Tyrosine Crystals
- Sulfa Crystals
- Indinavir Crystals

I. Acidic Urine Crystals

A. Amorphous Urates (Anhydrous uric acid)

- Normally present in urine in different quantity.
- Have pink to "brick red" color.
- From very small granules and seen in cluster.
- Dissolve in urine when the sample is gently heated.
- When urine is left in the refrigerator, it shows heavy precipitation of urates.



B. Uric Acid Crystals

- Polymorphs (different in shape) i.e. square, prism, hexagonal, etc.
- Yellow to yellow brown in color.
- Size is 30-150 μm
- Small quantity found in normal urine, but increases in association with:
- Increased Purine metabolism in case of gout.
- Increased Nucleic Acid turn over, such as leukemia.




Figure 3.8. Uric Acid Crystals

G. Bilirubin

- Very rarely seen.
- Have reddish brown color.
- Seen in case of elevated Bilirubin.
- Have various tiny squarish, beads or amorphous needle shape.
- Size is 5 μm (half RBC).
- Chemical test for bile pigments positive.

I. Acidic, Neutral, or Basic Urine Crystals

• Calcium Oxalate Crystal

- Are colorless and refractive.
- Have octahedral, envelope, shape.
- Size 10-12 μm .
- Normally seen in small amount.
- After consumption of high calcium, or oxalate rich foods, such as milk, tomatoes, asparagus, and orange, normally the crystals may be seen.



- In dehydration condition, such as, in hot weather where there is high perspiration and only a small amount of water is consumed per day Calcium oxalate crystals may be seen.
- Pathologically in large quantity may be seen in (severe chronic renal disease, and urinary calculus).




Figure 3.9. Calcium oxalate crystals

II. Alkaline, Neutral, or Slight Acidic Urine Crystals

• Triple Phosphates

- Colorless and refractive.
- Have “coffin lids” 3 to 4 to 6 – sided prism.
- Shape, or fern leaf or star shape.
- Size 130- 150 μ m.
- Seen in urine stasis (obstructive uropathy), or in urinary tract infections.
- Their presence is frequently indicative of bacterial infection by proteus mirabilis.






Fig.3.10. Triple phosphate crystals

III. Alkaline Urine Crystals

- ***Amorphous Phosphates***

- Normally seen in alkaline urine.
- Small, whitish granules usually seen scattered, & Soluble in 100g/1 acetic acid.

B. Calcium Carbonate

- Less commonly seen.
- Colorless.
- Have needle, spherical or dumbbells shape.
- Have very small crystals.
- If 100g/1, i.e.10% acetic acid is added, they dissolve, give off bubbles of gas.



Fig. 3.11. Calcium Carbonate crystals

C. Calcium Phosphates



- Seen in small amount in normal individual urine, and when they are in large amount, may indicate chronic cystitis, or prostatic hypertrophy.
- . Have star or needle shape.
- Colorless.




Fig. 3.12. Calcium Phosphates

MISCELLANEOUS

I. Spermatozoa

- Are small structures consisting of a head and tail, connected by a short middle piece (neck).
- Easily recognized especially if they are motile.
- Frequently seen in the urine of males.
- They may be seen in the urine of females, when the urine collected after coitus usually not reported, unless the physician has special interest in it.

II Mucus Trades

- Formed by the precipitation of mucoproteins in cooled urine.
- Normally little mucus trades seen in normal individuals.
- Have fine, fiber like appearance.
- Wavy in shape and tapered at ends.
- If not examined carefully may confuse with hyaline casts.
- Their presence in large amount with WBCs may indicate UTI.

III. Other Contaminates and Artifact Structure

- Muscle fibers
- Vegetable cells - all are fairly seen and easily
- Cotton fibers (wool fibers) recognizable.
- Structure from slide or cover slide - high retractile and non-uniform in size.

Fat droplets (other bubbles)

- Not evenly distributed.
 - Oil droplets
 - Pollen grains - are seasonal.
 - Starch granules - incomplete digestion of starch
- They can be confirmed by logon iodine.

*** To minimize the above mentioned contaminants and artifacts**

- Don't use dirty containers, slides and cover slides.



- Don't let urine specimen to open-air.
- Avoid contamination of urine with fats and oils.
- Avoid the drying of sediments.

Self-Check 5	Written Test
---------------------	---------------------

Instruction 1: Say True or False

1. The number of casts preserved decrease as the pH of the urine decreases.
2. Presence of RBCs in the urine is always indicative of a renal disease.
3. Waxy casts are the end stage in the degeneration of cellular casts.
4. Pyuria refers to elevated numbers of leucocytes in the urine.
5. The presence of Bacteria in the Urine is determined using only Microscope.

Instruction 2:- Choose the best possible answer for each of the following questions


1. The crystal found in the diagram below is-----?



- A. Triple Phosphate C. Amorphous Phosphate
 - B. Calcium oxalate D. Cysteine crystal
2. Which one is true about crystal in the Q-1 above?
 - A. Found in acidic urine only C. It is iatrogenic' crystals
 - B. Composed of magnesium, ammonium and phosphate D. All of the above
 3. The normal yellow color of urine is produced by:
 - A. Bilirubin B. Urochrome C. Urobilinogen D. Hemoglobin
 4. Ms. Darmi brought the urine specimen which is Yellow brown or "beer brown" in color. This may indicate?
 - A. presence of hemoglobin C. presence of bilirubin
 - B. Indicates hematuria D. Presence of protein
 5. Which of the following elements are commonly confused with RBCs?
 - A. Yeast cells B. leukocytes C. bubbles D. All
 6. In a microscopic examination of clear urine that produces a pink precipitate after refrigeration will show triple phosphate crystals.
 - A. True B. False C. Unknown



7. A yellow-brown specimen that produces a yellow foam when shaken can be suspected of containing:
A. Bilirubin B. Carrots C. Hemoglobin D. Rhubarb
8. The fungus *Candida* species (*Candida albicans*) are common in patients with
A. Diabetic mellitus C. Immunosuppressive therapy
B. Vaginites D. All of the above
9. What is the name of this crystal?



- A. Calcium phosphate C. Amorphous urate
B. Calcium oxalate D. Uric acid crystal
10. Which one is not true about crystal in Q-9 above?
A. Found in alkaline urine C. Indicate Presence of bilirubin in the urine
B. May found in monohydrate & dehydrate form D. All are true

Answer sheet

Instruction 1

1. _____
2. _____
3. _____
4. _____
5. _____

Instruction 2

1. _____
2. _____
3. _____
4. _____
5. _____
6. _____
7. _____
8. _____
9. _____
10. _____


Score _____
Rate _____



Remark _____

Name _____ ID.No _____ Date _____



 Operation sheet #4	LG #51	Performing Urine and Body Fluid analysis
	Topic	Performing Examination of Urine specimen

4.1. Performing urine examination

Procedure for Microscopic examination of Urine specimens

1. Collect all necessary materials used for the collection, centrifugation and examination of urine specimens
2. Preparation of patient
3. Explain the purpose of the test by using simple language. Do not use medical terms or try to explain details of the procedure.
4. Advise the patient how to collect the specimen. The first morning urine or mid-stream urine specimen is more preferable, because it is more concentrated.
5. If the patient is female, advice her to wash her genital organ before giving the specimen. This is because bacteria that are normally found on the genital tract may contaminate the sample and affect the result.
6. Advise the patient to collect at least 15 ml of urine in to the clean, sterilize and dry urine cup that is supplied from the laboratory.
7. The collected urine sample should arrive at a diagnostic laboratory as soon as possible.
8. Centrifugation of the urine specimen
9. Mix the urine specimen
10. Transfer about 10 ml of urine in the centrifuge tube.
11. Balance tubes in the centrifuge.
12. Centrifuge the specimen at a medium speed (from 1500 –2000 rpm) for 3-5 minutes
13. Discard the supernatant by quick inversion of the tube
14. Re suspend the sediment that is at the bottom of the tube, by tapping the tube by your fingers
15. Take the sediment by Pasteur pipette from the tube and transfer a drop into the clean, sterilized and dry slide. If Pasteur pipette is not available, gently incline the tube and place drop of sediment into the clean, sterilized and dry slide.
16. Apply cover slide on the urine sediment that is on the slide. This will make specimen to be spread on the slide on one cell thickness.
17. Put the slide on the stage of microscope and tie it by clips on the stage.
18. Lower the condenser, close the diaphragm and look under 10x objective of the microscope. Casts tend to concentrate near the edge of cover slide.
19. Then after looking through at least 20 fields of the low power objective, change the objective in to 40x objective. Do not forget to raise the condenser and opening of the diaphragm when you change the objective in to the high power (40x). Under high power objective also you should have to look for a minimum of 10-15 fields).



20. Then report what you get under 10 x (low power) and 40 x (high power) on the laboratory request form of the patient.

- *For determination of cellular elements, casts, etc, the number of elements seen under at least 10 fields should be counted and the average of this number is used for report value. Other elements such as parasites are usually reported as well.*


LAP TEST	PRACTICAL DEMONSTRATION
-----------------	--------------------------------

Name _____ ID.No _____ Date _____
Time started _____ Time ended _____

Instruction:- Demonstrate the following tasks(1hr)

1. Perform microscopic examination of urine according to the SOPs?
2. Identify urine crystals microscopically?



 Information sheet #6	LG #51	Performing Urine and Body Fluid analysis
	Topic	Performing Body Fluid Analysis

6.1. Performing Body Fluid Analysis

Objectives: At the end of this chapter the trainees be able to:

- ✓ Describe the overview of body fluids
- ✓ Describe body fluid analysis methods.
- ✓ Perform semen analysis.
- ✓ Perform cerebrospinal fluid analysis.

6.1. Cerebrospinal fluid (CSF)

Fluid in the space called sub-arachnoids' space between the arachnoids mater and pia mater

Protects the underlying tissues of the central nervous system (CNS)

- Serve as mechanical buffer to
- Prevent trauma,
- Regulate the volume of intracranial pressure
- Circulate nutrients
- Remove metabolic waste products from the CNS
- Act as lubricant

Has composition similar to plasma except that it has less protein, less glucose and more chloride ion

- It is one of the vertebrates body fluid contained in the cavity that surrounds the brain and the spinal cord.
- It supplies nutrients to the tissues of the central nervous system
- It helps to protect the brain and spinal cord from injury.
- The volume of the CSF in adults is 100–150ml; in children the volume is less and varies according to the body length.
- Maximum volume of CSF
 - Adults 150 mL
 - Neonates 60 mL
- Rate of formation in adult is 450-750 mL per day or 20 ml per hour
 - reabsorbed at the same rate to maintain constant volume
- Collection by lumbar puncture done by experienced medical personnel
- About 1-2ml of CSF is collected for examination
 - lumbar puncture is made from the space between the 4th and 5th



lumbar vertebrae under sterile conditions.

- Collected in three sequentially labeled tubes
 - Tube 1 Chemical and immunologic tests
 - Tube 2 Microbiology
 - Tube 3 Hematology (gross examination, total WBC & Diff)
 - This is the least likely to contain cells introduced by the puncture procedure



Fig. 3.13. location of CSF collection

Lab analysis

Clinical Significance

- Diagnosis of meningitis of bacterial, fungal, mycobacterial and amoebic origin or differential diagnosis of other infectious diseases
- subarachnoid hemorrhage or intracerebral hemorrhage

Principle of the test

- CSF specimen examined visually and microscopically and total number of cells can be counted and identified

Specimen: the third tube in the sequentially collected tubes*

- must be counted within 1 hour of collection (cells disintegrate rapidly). If delay is unavoidable store 2-8°C.
- All specimens should be handled as biologically hazardous

6.2. Semen analysis

- ❖ Used in the evaluation of reproductive dysfunction (infertility) in the male
- ❖ Used to select donors for therapeutic insemination



- ❖ Is a cost-effective and relatively simple procedure.
- ❖ Consists of microscopic and macroscopic components

Tests for semen

- ❖ Macroscopic
 - Physical (volume, viscosity, liquefaction)
 - chemical (eg. pH)
- ❖ Microscopic
 - stained preparation
 - wet-mount
- ❖ When investigating infertility, the basic analysis of semen (seminal fluid) usually includes:
 - Measurement of volume
 - Measurement of pH
 - Examination of a wet preparation to estimate the percentage of motile spermatozoa and viable forms and to look for cells and bacteria.
 - Sperm count
 - Examination of a stained preparation to estimate the percentage of spermatozoa with normal morphology.

Caution: Handle semen with care because it may

- contain infectious pathogens, e.g. HIV, hepatitis
- viruses, herpes viruses.

Macroscopic Examination

Measure the volume

- ❖ Normal semen is thick and viscous when ejaculated.
- ❖ It becomes liquefied usually within 60 minutes due to a fibrinolysin in the fluid.
- ❖ Failure to liquefy may indicate inadequate prostate secretion.
- ❖ When liquefied, measure the volume of fluid in millilitres using a small graduated cylinder.
- ❖ *Normal specimens:* Usually 2 ml or more

Measure the pH

- ❖ Using a narrow range pH paper, e.g. pH 6.4–8.0, spread a drop of liquefied semen on the paper.
- ❖ After 30 seconds, record the pH.
- ❖ *pH of normal semen:* Should be pH 7.2 – 7.8
- ❖ When the pH is over 7.8 this may be due to infection.



- ❖ When the pH is below 7.0 and the semen is found to contain no sperm, this may indicate dysgenesis (failure to develop) of the vas deferens, seminal vesicles, or epididymis.

Microscopic Examination

- ❖ be performed to obtain estimates of sperm concentration, motility, and agglutination.
- ❖ polygonal cells of the urethral tract and 'round cells' such as spermatogenic cells and leukocytes can also be observed when sperm are counted in a hemocytometer.
- ❖ Motility (normal range 50% or above) is expressed as the percentage of sperm that move.

- **Estimate the percentage of motile and viable spermatozoa**

Motility

- ❖ – Place 1 drop of *well-mixed* liquefied semen on a slide and cover with cover glass.
- ❖ – Focus the specimen using the 10_ objective.
- ❖ – Ensure the spermatozoa are evenly distributed
- ❖ – if not, re-mix the semen and examine anew preparation.
- ❖ – Using the 40_ objective, examine several fields
- to assess motility, i.e. whether excellent (rapid and progressive) or weak (slow and non progressive).
- Count a *Normal motility*: Over 50% of spermatozoa are motile within 60 minutes of ejaculation.

Reporting of results

- Motility (normal range 50% or above) is expressed as the percentage of sperm that move.
- Sperm moving rapidly in a straight line with little yaw and lateral movement are Grade 4
- if they move more slowly, Grade 3.
- Grade 2 sperm move even more slowly and with substantial yaw.
- Grade 1 sperm have no forward progression.
- Zero progression denotes absence of any motility
- If motility is less than 50%, a viability stain of eosin Y with nigrosin as a counterstain is done.
- dead sperm will stain red, whereas live sperm will exclude the dye and appear unstained.
- In samples with no visible sperm, such as post-vasectomy semen, the entire sample should be centrifuged, and the pellet examined for intact or



damaged sperm fragments.

- The spermatozoa remain motile for several hours.
- Perform gram stain smear:
 - When more than 60% of spermatozoa are non motile,
 - when more than a few leucocytes and
 - > 6 red blood cell/ HPF
 - Look for the type of bacteria that exist in the semen

Viability

procedure

- Mix one drop of semen with 1 drop of 0.5% eosin solution on a slide.
- After 2 minutes examine microscopically.
- Use the 40X objective to count the percentage of viable and non-viable spermatozoa.
- Viable spermatozoa remain unstained,
- non-viable spermatozoa stain red.
- *Normal viability: 75% or more of spermatozoa should be viable (unstained).*

Self-check 6	Written examination
--------------	---------------------

Instruction1:-Say true or false for each of the following questions


1. Semen analysis used in the evaluation of reproductive dysfunction (infertility) in the male
2. CSF is used to select donors for therapeutic insemination
3. Lumbar puncture is made from the space between the 4th and 5th lumbar vertebrae under sterile conditions.
4. CSF should be collected in three sequentially labeled tubes among them Tube 1 is used for Chemical and immunologic tests.
5. Tube 2 is used for Hematology (gross examination, total WBC & Diff) whereas Tube 3 is for Microbiology tests.

Answer sheet

1. _____
2. _____
3. _____
4. _____
5. _____

Score _____
Remark _____



 Operation sheet #5	LG #51	Performing Urine and Body Fluid analysis
	Topic	Procedures for Collection and transportation of semen

6.1. Procedures for Collection and transportation of semen

1. Give the person a clean, dry, leak-proof container,
and request him to collect a specimen of semen at home following 3 days of sexual abstinence. Condom is used to collect the fluid, this must be well-washed to remove the powder which coats the rubber. It must be dried completely before being used.
2. Label the container (name, date and time of collection, period of abstinence)
3. Deliver the specimen to the laboratory within 1 hour
4. Fluid should be kept as near as possible to body temperature.
5. This is best achieved by placing the container inside a plastic bag and transporting it in the person's armpit . .
 - a. **Procedure for Estimating the percentage of motility of spermatozoa**
 1. Place 1 drop of *well-mixed* liquefied semen on a slide and cover with cover glass.
 2. Focus the specimen using the 10_ objective.
 3. Ensure the spermatozoa are evenly distributed
 4. if not, re-mix the semen and examine anew preparation.
 5. Using the 40_ objective, examine several fields
 - b. **Procedure for Estimating the percentage of viability of spermatozoa**
 1. Mix one drop of semen with 1 drop of 0.5% eosin solution on a slide.
 2. After 2 minutes examine microscopically.
 3. Use the 40X objective to count the percentage of viable and non-viable spermatozoa.
 4. Viable spermatozoa remain unstained,
 5. non-viable spermatozoa stain red.
 6. *Normal viability: 75% or more of spermatozoa should be viable (unstained).*



Lap test	Practical demonstration
-----------------	--------------------------------

Name _____ ID.No _____ Date _____
Time started _____ Time ended _____

Instruction:- Demonstrate the following tasks(1hr)


Project 1:- Performing semen analysis

Task1:- Perform Collection and transportation of semen

Task2:- Estimate the percentage of motility of spermatozoa

Task3:- Estimate the percentage of viability of spermatozoa



 Information sheet #7	LG #51	Performing Urine and Body Fluid analysis
	Topic	Quality control in urinalysis

Quality control in urinalysis.

- ❖ Quality assurance is a set of activities starting from specimen collection to issuing test results that ensure test results are accurate and precise as possible.
- ❖ It is the sum of all the activities of the laboratory that ensures test results are of good quality.
- ❖ Quality assurance includes
 - ❖ inside and outside the laboratory performance standards
 - ❖ good laboratory practice and management skills that are required by achieving and maintaining a quality service and that provide for continuing improvement
 - Part of quality assurance, which primarily concern the control of errors in the performance of tests and verification of test results.
 - must be practical, achievable, affordable, and above all continuous
 - The purpose of quality control procedure is to monitor analytical processes, analytical error and to correct result of analysis.

Two types of quality control programs

A) Internal quality control

- Is carried out in the laboratory, an intra-lab program.
- Encompasses all measurements made, technical skills performed within an individual laboratory.
- use control samples, like pooled serum
- The purpose of quality control program is to insure tests are performed reliably and reported correctly.
- Effective quality control systems detect errors at an early stage, before they lead to incorrect test results.

■ B) External quality control.

1. External quality control is observation of variance in results when the same material is analyzed in different laboratories

External quality control is observation of variance in results when the same material is analyzed in different laboratories

Quality control steps:

- Pre analytical steps



- Analytical steps

- Post analytical steps

1. Pre analytical Quality control in urinalysis

- Read and understand requested paper
- guide the patient to bring an appropriate urine sample
- Labeling the urine container after collecting the sample
- Check the material we are going to use whether they are properly cleaned or not
- Ask the patient whether the urine sample left ,more than two hours, after it is voided.
- Do not accept contaminated requested paper
- Check the slide, the microscope, & all needed material before taking the next procedure.
- If the urine comes from far place ask or read the preservative applied
- Concentrate and find out an abnormality related to chemical & physical appearance.
- Proper sample preparation is also most important.
- Reduce possible source of errors
- Do not open the centrifuge while it is not stopped
- Proper balance of urine in the centrifuge

2. analytical quality control in urinalysis

- small urine sample how to be rejected
- Follow exactly standard operation procedure (sop)
- Check and read reagent strip chemical test according to the instruction of the manual of the manufacturer, at the right time
- Write the physical appearances properly
- use the needed amount of urine for centrifugation
- When discarding the supernatant, it has to be quick and vertical upside down in order not to lose the sediment
- Examine as quickly as possible

3. post analytical quality control in urinalysis

- Proper written result
- Correct calculation
- Result interpretation



Self-Check 7	Written Test
---------------------	---------------------

Exercises4: Say True or False

1. External quality control is observation of variance in results when the same material is analyzed in different laboratories
2. Quality assurance is a set of activates starting from specimen collection to issuing test results that ensure test results are accurate and precise as possible.
3. Quality control is the sum of all the activates of the laboratory that ensures test results are of good quality.
4. Quality assurance includes inside and outside the laboratory performance standards
5. Proper sample preparation is part of post-analytical quality assurance.

Answer sheet

Name _____ ID.No _____

1. _____
2. _____
3. _____
4. _____
5. _____

Score _____
Remark _____



3.1. Verifying Laboratory results before releasing for clinician/client

In this topic, a review is given of all elements involved in the validation of clinical laboratory results. Validation will include:

1. Method validation,
2. Instrument validation,
3. Preanalytical validation procedures,
4. Analytical validation procedures, and
5. Postanalytical validation procedures.

Within the scope of this sub-topic, all of these different elements are discussed in detail. The management of all these steps is the only way to guarantee a correct result, if this is used either for patient treatment or in clinical evaluation studies.

All the types of validation is expressed in the diagram in page 64 below

Checklist for validation of test results

A validation of patient results should be performed using this checklist. Only when a complete validation is performed the report may be authorized to be sent to the requester.

Patient ID: _____

Pre-analytical phase

- ✓ Patient was correctly identified
- ✓ Patient was properly prepared for sample collection
- ✓ The person collecting the samples was correctly identified
- ✓ Sample was labeled correctly and clear
- ✓ The request form matches the specimen
- ✓ The request form contains correct and clear contact details of the requester
- ✓ The date and time of collection is indicated on the request form
- ✓ The specimen was transported appropriately to the laboratory
- ✓ The specimen was received in acceptable condition
- ✓ The log book entry matches the specimen label

Analytical phase

- ✓ Reagents and test kits used were within expiry date
- ✓ Quality controls associated with the result were acceptable
- ✓ There were no flags on the analyzer's results that need investigation



- ✓ If diluted, the final results were calculated correctly with the correct dilution factor
- ✓ Results are within the biological reference intervals
- ✓ Panic (critical) values are confirmed
- ✓ The results make clinical sense
- ✓ Confirmatory testing or established testing algorithms were completed
- ✓ If applicable: previous patient results are available to assist with interpretation of current sample's result

Post Analytical phase

- ✓ The report shows an appropriate result including test and result match for each test requested
- ✓ Proper concentration units for results are used
- ✓ The decimal place is correct (if results have decimals)
- ✓ The persons performing the tests are identified
- ✓ All results and documentation are legible
- ✓ In case of results within critical intervals the need for immediate notification is indicated on the report and an immediate notification form is used to verify correct reception of the result report by the requester
- ✓ If applicable, the report contains interpretative information to assist the clinician
- ✓ The release of the results is dated and timed


Remarks: _____

Authorizer's name, signature and date for completion of validation and correctness of results:

Date: _____ Name: _____ Signature: _____



validation of results





References

1. Cheesbrough, M.(2005)Medical Laboratory Manual for Tropical Countries. Vol.1 Cambridge: Butterworth-Heinemann. Ltd
2. Free A.H (1975) Urinalysis in Clinical Laboratory Practice. Cleveland: CRC Press
3. Sood, Ramnik(1987)Medical Laboratory Technology Methods and Interpretation. 4thed. New Delhi India: Jaypee Brothers
4. Linne, J.J and Ringsru, K.M (1979). Basic technique for the medical laboratory, McGraw-Hill,inc
5. Assamenew K., Mistire W., Belayhun K. (2002) Urinalysis for Medical Laboratory Technology Students.
6. Brunzel, N.A. (1994). Fundamentals of urine and body fluid analysis. Philadelphia: W.B. Saunders

Prepared By							
No	Name	Educational Background	LEVEL	Region	College	Email	Phone Number
1	Kalicha Boru	Laboratory	B	oromia	Nagelle HSC	boru9683@gmail.com	0912493885
2	Furo Beshir	Laboratory	A	Harari	Harar HSC	nebi.furo@gmail.com	0911739970
3	Motuma Chali	Laboratory	B	oromia	Nekemte HSC	lammiifcaalii@gmail.com	0938456753
4	Abdirahman Mahad	Laboratory	A	Somali	Jigjiga HSC	abdirahman7584@gmail.com	0911044715
5	Adisu Tesfaye	Laboratory	B	Somali	Jigjiga HSC	adistesfaye21@gmail.com	0931747320
6	Kebebe Tadesse	Laboratory	B	BGRS	Pawi HSC	no	0926841290
7	Tagel Getachew	Laboratory	A	Harari	Harar HSC	tagegetachew@gmail.com	0915746748