



Medical Laboratory

NTQF Level III

Learning Guide 22

Unit of Competence: Collect and Process Medical Samples

Module Title: Collect and Process Medical Samples

LG Code: HLT MLS4_LO1-08

TTLM Code: HLT MLS4TTLM 0919v1

LO1. Apply concept of physiology and anatomy |



Learning Instructions:

1. Read the specific objectives of this Learning Guide.
2. Follow the instructions described in number 5 to 11.
3. Read the information written in all Information Sheets. Try to understand what are being discussed. Ask you teacher for assistance if you have hard time understanding them.
4. Accomplish all Self-check according to learning session separately
5. Ask from your teacher the key to correction (key answers) or you can request your teacher to correct your work. (You are to get the key answer only after you finished answering the Self-checks).
6. If you earned a satisfactory evaluation proceed to next Information Sheet. However, if your rating is unsatisfactory, see your teacher for further instructions or go back to pervious Learning Activity.
7. Submit your accomplished Self-check. This will form part of your training portfolio.
8. Try to understand what are being discussed. Ask you teacher for assistance if you have hard time understanding them.
9. Accomplish the all Self-checks.
10. After you accomplish self check proceed to operation sheet if available.
11. If you perform operation procedure correctly proceeds to LAP. test if available



Information Sheet-1	Apply concept of physiology and anatomy
---------------------	---

Learning outcomes (objectives):

At the end of this module the trainee will be able to:-

- Identify human anatomy and physiology
- Understand the type and nature of samples
- Identify the time of sample collection and collection sites

1.1. Concepts of human anatomy and physiology

1.1.1. What are Anatomy and Physiology?

- **Anatomy:** the word anatomy is derived from a Greek word “Anatome” meaning to cut up. It is the study of structures that make up the body and how those structures relate with each other.

The study of anatomy includes many sub specialties. These are Gross anatomy, Microscopic anatomy, Developmental anatomy and Embryology.

Gross anatomy studies body structure without microscope. Systemic anatomy studies functional relationships of organs within a system whereas Regional anatomy studies body part regionally. Both systemic and regional approaches may be used to study gross anatomy.

Microscopic anatomy (Histology) requires the use of microscope to study tissues that form the various organs of the body.

- ✓ **Gross (macroscopic) anatomy:** the study of structures large enough to be seen with the naked eye
- ✓ **Regional anatomy:** all the body structures (muscles, bones, blood vessels, nerves, etc.) in a given body region, such as the abdomen or leg, are examined at the same time
- ✓ **Systemic anatomy:** body is studied system by system
 - ✓ Example: when studying the cardiovascular system, you would examine the heart and the blood vessels of the entire body
- ✓ **Surface anatomy:** internal body structures as they relate to the overlying skin
 - ✓ Used when identifying the bulging muscles beneath a bodybuilder’s skin, and clinicians use it to locate appropriate blood vessels in which to feel pulses and draw blood
- ✓ **Microscopic anatomy:** the study of structures that are too small to be seen with the naked eye
 - Cytology: study of individual cells



-Histology: study of tissues

- ✓ Developmental anatomy: the study of the change in body structures over the course of a lifetime
 - ✓ Embryology: concerns developmental changes that occur before birth
- Physiology: the word physiology derived from a Greek word for study of nature. It is the study of how the body and its part work or function.
 - ✓ Renal physiology: concerns kidney function and urine production
 - ✓ Neurophysiology: explains the workings of the nervous system
 - ✓ Cardiovascular physiology: examines the operation of the heart and blood vessels

While anatomy provides us with a static image of the body's architecture, physiology reveals the body's dynamic nature

1.1.1. Levels of Structural Organization


- Chemical level is the simplest level of organization:
 - Atoms, tiny building blocks of matter, combine to form molecules such as water and proteins
 - Molecules combine in specific ways to form organelles, which are the basic unit of living cells
 - Cells are the smallest units of living things
 - All cells have some common functions, but individual cells vary widely in size and shape, reflecting their unique functions in the body.
- Cellular level: smallest unit of life, and varies widely in size and shape according to the cell's function
- Tissue level: groups of similar cells having a common function
 - Four basic tissue types: each tissue type has a characteristic role in the body
 - Epithelium: covers the body surface and lines its cavities
 - Muscle: provides movement
 - Connective: supports and protects body organs
 - Nervous: provides a means of rapid internal communication by transmitting electrical impulses
- Organ level: made up of discrete structures that are composed of a least two groups of tissues that work together to perform a specific function in the body
 - Stomach: epithelium lining, muscles, blood vessels, connective tissues, nerve fibers, etc.
- Organ system level: a group of organs that work closely together to accomplish a specific purpose
 - Respiratory and circulatory system, digestive and circulatory systems
- Organismal level: the total of all structures working together to promote life
 - The living human being



1.2. Maintaining Life Necessary Life Functions

- Maintaining Boundaries: allows an organism to maintain separate internal and external environments, or separate internal chemical environments
 - Integumentary System or Skin
- Movement: allows the organism to travel through the environment, and allows transport of molecules within the organism
 - Skeletal, Circulatory, Muscular Systems
- Responsiveness: or irritability, is the ability to detect changes in the internal or external environment and respond to them
 - Muscular System


ORGAN SYSTEMS



(a) Integumentary System
Forms the external body covering; protects deeper tissues from injury; synthesizes vitamin D; site of cutaneous (pain, pressure, etc.) receptors, and sweat and oil glands.



(b) Skeletal System
Protects and supports body organs; provides a framework the muscles use to cause movement; blood cells are formed within bones; stores minerals.



(c) Muscular System
Allows manipulation of the environment, locomotion, and facial expression; maintains posture; produces heat.

Copyright © 2004 Pearson Education, Inc., publishing as Benjamin Cummings.


Fg. 1.1. Integumenatry

Fg. 1.2.skeletal

Fg. 1.3. Muscularsystem


- Nervous System:
 - Responsiveness to external and internal environments by activating muscles and glands
- Endocrine System:
 - Regulating body functions such as: growth, reproduction, and nutrition
- Cardiovascular System:
 - Transportation of nutrients, waste, gases, and hormones throughout the body

ORGAN SYSTEMS



(d) Nervous System
Fast-acting control system of the body; responds to internal and external changes by activating appropriate muscles and glands.

Copyright © 2004 Pearson Education, Inc., publishing as Benjamin Cummings.



(e) Endocrine System
Glands secrete hormones that regulate processes such as growth, reproduction and nutrient use (metabolism) by body cells.



(f) Cardiovascular System
Blood vessels transport blood, which carries oxygen, carbon dioxide, nutrients, wastes, etc.; the heart pumps blood.


Fig. 1.4. Nervous

Fig. 1.5. Endocrine

Fig. 1.6. Cardiovascular system


- (Lymphatic System/Immunity:
 - Body defenses
- Respiratory System:
 - External and internal gas exchanges
- Digestive System:
 - Breakdown and absorption of nutrients

ORGAN SYSTEMS



(g) Lymphatic System/Immunity
Picks up fluid leaked from blood vessels and returns it to blood; disposes of debris in the lymphatic stream; houses white blood cells (lymphocytes) involved in immunity. The immune response mounts the attack against foreign substances within the body.

Copyright © 2004 Pearson Education, Inc., publishing as Benjamin Cummings.



(h) Respiratory System
Keeps blood constantly supplied with oxygen and removes carbon dioxide; the gaseous exchanges occur through the walls of the air sacs of the lungs.



(i) Digestive System
Breaks down food into absorbable units that enter the blood for distribution to body cells; indigestible foodstuffs are eliminated as feces.

Fig. 1.7. (Lymphatic


Fig. 1.8. Respiratory

Fig. 1.9. Digestive system




- Urinary System:
 - Absorption of waste from the blood and elimination
- (k): Male Reproductive System:
 - Production of sperm
- (l): Female reproductive System:
 - Production of eggs

ORGAN SYSTEMS




(j) Urinary System

Eliminates nitrogenous wastes from the body; regulates water, electrolyte and acid-base balance of the blood.



(k) Male Reproductive System

Overall function is production of offspring. Testes produce sperm and male sex hormone; ducts and glands aid in delivery of sperm to the female reproductive tract. Ovaries produce eggs and female sex hormones; remaining structures serve as sites for fertilization and development of the fetus. Mammary glands of female breasts produce milk to nourish the newborn.



(l) Female Reproductive System

Copyright © 2004 Pearson Education, Inc., publishing as Benjamin Cummings.

Fig. 1.10. Urinary

Fig. 1.11. Male reproductive

Fig. 1.12. Female reproductive system



Self-Check 1	Written Test
--------------	--------------

Answer the Following Questions (2 point each):

1. _____ is the science that study of structures that make up the body and how those structures relate with each other.

- | | |
|---------------|--------------|
| A. Anatomy | C. Biology |
| B. Physiology | D. Chemistry |

2. The System allows the organism to travel through the environment, and allows transport of molecules within the organism

- | | |
|--------------|---------------------------------|
| A. Urinary | D. Integumentary System or Skin |
| B. Digestive | |
| C. Nervous | |

3. Groups of similar cells having a common function

- | | |
|----------|-------------|
| A. Cell | C. Tissue |
| B. Organ | D. Organism |

4. Types of science used for study of nature. It is the study of how the body and its part work or function

- | | |
|---------------|-----------|
| A. Anatomy | C.A and B |
| B. Physiology | D. None |

5. List four basic tissue types

- | | |
|----------|----------|
| A. _____ | C. _____ |
| B. _____ | D. _____ |

Note: Satisfactory rating - 10 points Unsatisfactory - below 10 points

You can ask you teacher for the copy of the correct answers

Answer Sheet

Medical laboratory Level III	Version.01 Sep-2019	Page 8 of 20
	Author –Federal TVET Agency	



Score _____
Rating _____

Name: _____ Date: _____

Short Answer Question 1. _____ 2. _____ 3. _____ 4. _____

5. _____, _____, _____



Information Sheet-2	Types of specimen and purpose of sample
---------------------	---

1.2. Type and nature of samples are identified

1.2.1. Types of specimen

There are different types of common clinical specimen collected and/or analyzed in the diagnostic laboratory. These should include:

- Stool
- Urine
- Body discharges
- Blood
- Skin slip, skin slit and skin scrapping
- Cerebrospinal Fluid(CSF)
- Swabs (eg .throat swab)



COLLECTING SAMPLES FOR CLINICAL TESTING	SAMPLE TYPES
Routine clinical samples	Whole Blood Plasma Serum Red Blood Cells (RBC) Urine
Non-Routine Liquid Biological Samples	CSF Breast Milk Exudates Lavage Fluid Oral Fluid
Non-Routine Solid Biological Samples	Stool Meconium Hair Nails
Samples - Metals Analyses	Whole Blood Plasma Red Blood Cells (RBC) Serum Urine

Table 1.1. Types of specimen

Notes about Shipping Temperature: Ship specimens chilled unless indicated otherwise in the Labs Test procedure SOP.

Notes About Specimen Labeling: Please use water-resistant ink when labeling specimens. Include the following information on each specimen:

1. Name of the patient (subject, employee or patient)
2. ID# or case number
3. Specimen type (blood, serum, urine, hair, etc.)
4. Date and time specimen was obtained

Collecting Samples for Clinical Testing

The following are general considerations for specimen collection, appropriate collection containers for specific testing and guidelines for sample submission.



For trace element, volatile compound analysis, specimen types other than blood, serum, plasma or urine, read additional instructions.

Collecting Routine Clinical Specimens

Note: The use of serum/plasma separator tubes is not recommended for use when collecting specimens for testing to be performed at NMS Labs. The use of these tubes can compromise test results.

Whole Blood refers to whole blood drawn into evacuated tubes (e.g., Vacutainer® tube) containing anticoagulant.

Plasma should be obtained by drawing blood into an anticoagulant, evacuated tube. Invert tube at least eight times to mix, immediately centrifuge for approximately 10 minutes at 3000 rpm, and carefully transfer the supernatant plasma into a labeled plastic container (polyethylene vial).

Serum should be obtained by drawing blood into an evacuated tube containing no anticoagulant and with a non-coated interior. Allow 20 minutes for clotting. Immediately centrifuge for approximately 10 minutes at 3000 rpm and carefully transfer the supernatant serum into a labeled, plastic container (polyethylene vial).

- Red Blood Cells (RBC) should be prepared from whole blood using the procedure for plasma specimens. Transfer the plasma from the centrifuged red blood cell fraction and submit the red blood cells, labeling the specimen accordingly.
- Urine specimens should be collected in a plastic (polyethylene) bottle. Do not fill bottle past the shoulder. Close with a screw cap, which has a self-sealing liner (e.g., 50ml bottle and cap).
- Collecting Non-Routine Liquid Biological Specimens

While not all biological fluids are found in large volumes, make every effort to collect as much as possible. Place fluid into a labeled screw-capped plastic container for shipping.

- Collecting Non-Routine Solid Biological Specimens
 - ✓ Stool: Collect an entire bowel movement and place into a clean, labeled 40 mL polyethylene bottle. Fill no more than 2/3 of the bottle. Write the total original weight on the label (mass or volume) of the specimen. Freeze. View frozen sample shipping instructions.

Medical laboratory Level III	Version.01 Sep-2019	Page 12 of 20
	Author –Federal TVET Agency	



- ✓ Meconium: Collect at least 5 grams, approximately 1 tablespoon, of the black-tarry Meconium sample and place into a clean 40 ml polyethylene bottle. The sample may be combined several times from each evacuation up to approximately 72 hours or when the sample starts to turn yellowish-green.

- **Collecting Samples for Metals Analyses**

Specimens collected for trace-metals analyses must be protected from contamination during collection. Conduct specimen collection in a clean, dust free environment using appropriate certified metal-free collection containers. To reduce specimen contamination, powder-free gloves are recommended during collection. It is recommended that unexpected elevated results be verified by testing another specimen.

Blood, Plasma, or Red Blood Cells (RBC): To flush the collection device of metal contaminants, always draw a plain red top tube (no additive) first. Secondly, draw 7 mL of blood into a trace-metal free, EDTA Royal Blue Top Tube, unless otherwise indicated in the online test catalog specimen collection instructions. Always check the online test catalog for the most up to date collection requirements. Invert the Blue Top tube at least 8 times to ensure adequate mixing of the EDTA anticoagulant to prevent clotting. Discard the first red top tube appropriately unless being used for non-metals testing.

For RBC or Plasma samples, centrifuge the blood sample as soon as possible after collection for approximately 10 minutes at 3000 rpm. Place plasma into an acid washed plastic screw capped vial. Leave the RBCs in the original tube and replace stopper.

- ✓ Serum: To flush the collection device of metal contaminants, always draw a plain Red Top Tube (no additive) first. Secondly, draw 7 ml of blood into a trace-metal free, evacuated tube containing no anticoagulant. Refer to the online test catalog for the most up to date collection requirements. Allow the specimen to clot for at least twenty (20 minutes). Centrifuge for approximately 10 minutes at 3000 rpm and transfer all of the serum into an acid washed plastic screw capped vial. Discard the first red top tube appropriately unless being used for non-metals testing.
- ✓ Urine: Collect 20 mL of urine directly into a labeled trace-metal free (acid-washed, deionized water-rinsed, air-dried), polyethylene bottle with a plastic lined screw cap or an acid washed collection container. Always check the online test catalog for the most up to date collection requirements.
- ✓ Nails: Nail clippings may contain analytes of interest, which were deposited during the growth of the nail. Nail clippings (at least 0.5 g) from each finger or toe should be collected and packaged separately in acid washed plastic bottles. Each bottle should

Medical laboratory Level III	Version.01 Sep-2019	Page 13 of 20
	Author –Federal TVET Agency	



be labeled with the weight, if known, of the nail collected and its source, e.g., right index finger.



Self-Check 2	Written Test
--------------	--------------

Answer the Following Questions (2 point each):

1. List common clinical specimen used in laboratory
2. Types of blood specimens obtained by drawing blood into an anticoagulant _____
3. Specimen labeling include the information

Note: Satisfactory rating - 6 points Unsatisfactory - below 6 points

You can ask you teacher for the copy of the correct answers

Answer Sheet

Score _____
Rating _____

Name: _____

Date: _____

Short Answer Question

1. _____

2. _____

3. _____



Information Sheet-3	Time of sample collection and collection sites
---------------------	--

1.3. Time of sample collection and collection site

1.3.1. Time of collection:

The time of collection provide best chance of recovery of the causative agent .For sputum and urine the preferred time is early in the morning soon after the patients awaken.

For blood specimen the time of collection should be when the patient's temperature begins to rise.

- Collection of specimens before the administration of antimicrobial:

Because antimicrobials limits recovery of pathogens if specimen is collected after the administration of anti-microbial.

- Age of specimens:

Age of the specimen directly influences the recovery of protozoan organism

- Stage of the disease at which the specimen is collected

Enteric pathogens are present in great numbers during the acute or diarrheal stage of intestinal infection

1.3.2. Site selection consideration:

- Clinician should locate right anatomic site & select appropriate tests &specimens based on:
 - ✓ Physical examination (sign & symptoms)
 - ✓ Radiological examination
- Laboratory personnel should collect specimens from actual infection site with little external contamination by using aseptic technique and sterile container and should collect specimens from right site:
 - ✓ To prevent contamination of specimen &
 - ✓ To protect the patient from infection
- Sites of Infection where the Specimen is likely to become contaminated during collection
 - ✓ Sample from lower respiratory tract can be contaminated from Oro-pharynx
 - ✓ Sample from bladder can be contaminated from urethra
 - ✓ Sample from cervix can be contaminated from vagina

Medical laboratory Level III	Version.01 Sep-2019	Page 16 of 20
	Author –Federal TVET Agency	







Self-Check 3	Written Test
--------------	--------------

Answer the Following Questions (2 point each):

1. For sputum and urine the preferred time is specimen collection is
 - A. Early in the morning soon after the patients awoken.
 - B. At any time of the day
 - C. Random
2. For blood specimen the time of collection should be when the patient's.
 - A. Only at Early morning
 - B. After treatment
 - C. When body temperature begins to rise
 - D. At any time
3. Why laboratory personnel should collect specimens from actual infection site?

Note: Satisfactory rating - 6 points Unsatisfactory - below 6 points
You can ask you teacher for the copy of the correct answers

Answer Sheet

Score _____
Rating _____

Name: _____

Date: _____

1. _____
 2. _____
 3. _____
- _____
- _____

Prepared By

No	Name	Educational Back grund	LEVEL	Region	College	Email	Phone Number
1	Kalicha Boru	Laboratory	B	oromia	Nagelle HSC	boru9683@gmail.com	0912493885
2	Furo Beshir	Laboratory	A	Harari	Harar HSC	nebi.furo@gmail.com	0911739970
3	Motuma Chali	Laboratory	B	oromia	Nekemte HSC	lammiifcaalii@gmail.com	0938456753
4	Abdirahman Mahad	Laboratory	A	Somali	Jigjiga HSC	abdirahman7584@gmail.com	0911044715
5	Adisu Tesfaye	Laboratory	B	Somali	Jigjiga HSC	adistesfaye21@gmail.com	0931747320
6	Kebebe Tadesse	Laboratory	B	BGRS	Pawi HSC	kebebetadesse@gmail.com	0926841290
7	Tagel Getachew	Laboratory	A	Harari	Harar HSC	tagegetachew@gmail.com	0915746748